Mini L.E.D.

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ENGLISH

SUMMARY

1- SAFETY INSTRUCTIONS	4
2- PRESENTATION	5
3- INSTALLATION	5
4- PRODUCT DESCRIPTION	5
5- MODES USED	6
6- EVERYDAY USE	6
7- CARE	6
8- TROUBLESHOOTING	7
9- SPECIFICATIONS	8
10- ELECTROMAGNETIC COMPATIBILITY	9
11- DISPOSAL AND RECYCLING	11
12- LIABILITY	11
13- REGULATIONS	11
14- SYMBOLS AND ABBREVIATIONS	12

1 - SAFETY INSTRUCTIONS

Warning: Before using your **MINI L.E.D.**, please read the whole of this Safety Instructions section carefully. The

manufacturer's warranty is applicable only if these instructions concerning the operation and the safety of the unit have been properly applied. All these safety measures assume a good level of knowledge of dental techniques and photo-polymerization, and the more specific knowledge regarding the use of the **MINI L.E.D.** contained in this Operating Manual.

- 1) This unit must be used in strict accordance with the instructions contained in this Operating Manual. The manufacturer declines any liability if these rules are not followed or if the unit is used for any other application.
- 2) Before connecting the unit, make sure that the power voltage is compatible with the voltage marked on the power adapter used to charge the curing light. A different voltage would damage the unit and could result in injury to the patient and/or the operator.
- 3) The light radiation produced by this type of unit can be dangerous and must not be pointed at the eyes, even if the practitioner or the patient is wearing suitable protective eyewear (class 2M laser). The light produced by this equipment must be directed only at the zone to be treated in the oral cavity.
- 4) Any condensation inside electrical equipment can be a hazard. If the curing light has to be moved from a cool place to a warm place, do not use it immediately: let it reach the new room temperature first.
- 5) Do not replace the battery during use.
- 6) Comply with all the instructions given in this Manual. Do not insert or attempt to insert metal objects into the unit or its base station, as this could cause an electrical discharge.
- 7) The manufacturer declines all liability if damaged parts or accessories are not replaced exclusively by the manufacturer's products. In particular, the use of other light guides, power adapters or batteries could be dangerous for the patient and the operator.
- 8) In the case of a fault, unplug the base station immediately and make sure that no one can use the curing light before it is checked by the manufacturer or dealer. The fault may be the result of non-compliance with the safety rules or of a technical problem in the unit.
- 9) Do not use the unit near a heat source. Do not use solvents, detergents or flammable products, which may damage the unit or cause short-circuits.
- 10) Only the unit manufacturer or dealer is qualified to carry out repairs.
- 11) All lights used for photo-polymerization, including the **MINI L.E.D.**, must not be used on persons who are suffering or have suffered from photo-biological reactions (including urticaria solaris or erythropoietic protoporphyria) or who are receiving treatment including photosensitizing medications (including methoxsalens or chlorotetracycline).
- 12) Persons (practitioners or patients) who have previously suffered from a retina or crystalline lens condition or who have undergone eye surgery, in particular for cataracts, must see their ophthalmologist before using the **MINI L.E.D.**. Even if the eye specialist agrees, prudence is strongly recommended, because the intensity of the light could cause accidents. Class 2M protective eyewear, suitable for use with units emitting radiation at wavelengths between 420 and 480 nm, should be worn at all times.
- 13) **MINI L.E.D.** must not be used if the patient and/or the operator has a cardiac stimulator or any other active implant (e.g. a cochlear implant).
- 14) The device is not designed to withstand shocks delivered by an electric defibrillator.
- 15) Any variation of the power supply voltage or the magnetic field outside the limits in force could switch the curing light to automatic operation or interfere with its operation.
- 16) Electromagnetic interference: the unit complies with the applicable standards (IEC 60 601-1-2) for emission (radiated electrical field and interference with the supply voltage) and for immunity (protection against electrical fields, fast supply transients in bursts, electrostatic discharge and supply surges).
- 17) Your dealer carries stocks of all spare parts that you might need. Used batteries must be returned to the dealer.

- 18) For use only by dental professionals.
- 19) When you need to transport the unit, unscrew the battery and protect the light guide against impacts.
- 20) For long-term storage we recommend that you remove the battery from the unit body so that it does not discharge slowly and risk damage.
- 21) When your unit reaches the end of its life, we recommend that you contact your dental equipment dealer or, failing this, one of the ACTEON sites (listed at the end of this manual) to find out what procedure is advised.

2 - PRESENTATION

The **MINI L.E.D.** produces a visible blue light between a wavelength of 420 and 480nm for the photo-polymerization of dental materials. After checking the packaging for damage (please retain the packaging), you will find the following:

- the handpiece, fitted with a protective cap, incorporating the light-emitting diode,
- the base station for recharging the unit, with a built-in radiometer,
- the power adapter,
- the lithium ion battery,
- the multifibre light guide, 45° curvature, Ø7.5mm, sterilizable,
- the protective light shield,
- the Operating Manual and other documents.

Available options:

- a multifibre light guide, opalescent, curvature 45°, Ø5.5mm, sterilizable

3 - INSTALLATION


- 1- Prior to each use, sterilize the light guide and disinfect both base station and handpiece (see also the chapter on care).
- 2- Screw the battery onto the handpiece, insert the sterilized light guide into the handpiece. Ensure that the light guide is inserted correctly (the friction contact makes a click).
- 3- Place the base station on a stable surface. Plug in the mains adapter and ensure that the plug is suitable for the voltage of your power supply. The green light will

flash twice.

- 4- Place the handpiece on its base station (there will be an audible beep each time the lamp is put in the base) and leave it to charge for at least 2 hours 30 minutes, to make sure that the battery is fully charged.
- 5- The green light starts to flash indicating that the battery is being charged. As soon as the battery is fully charged, the light stops flashing and stays lit.

4 - PRODUCT DESCRIPTION

Now that the **MINI L.E.D.** is fully charged, you need to familiarize yourself with its buttons and menus:

Handpiece (1)

a) LED indicators:

The handpiece has four LED indicators:

- One LED located at the top, next to the light guide, indicates the status:

Color	Mode
Green	Normal mode
Red	Power supply problem/battery low
Red flashing	Thermal protection

- A row of three LEDs indicates the curing mode selected.

Symbol	Mode
	Fast
	Pulse
	Stepped

b) Buttons

The handpiece has two buttons:

- one navigation button which is used to switch from one

menu to another

- one on/off light button (trigger) close to the light guide.

Base station (2)

a) Wattmeter

A built-in wattmeter, associated with an LED indicator allows the operator to check accurately the curing light's performance. The aperture is located on the front part of the base station. An adjustable ring allows you to adapt the wattmeter's input aperture to the Ø5.5mm light guide.

b) To proceed

- Set the device in Fast Curing mode.
- Insert the light guide into the handpiece after checking no VLC material remains on the light guide.
- Place the end of the light guide in the aperture and activate the **MINI L.E.D.**.

c) LED indicator

Color	Mode
Green	Power over 800mW/cm ²
Red	See troubleshooting chapter

d) Charge

When the handpiece is inserted in the base station, another LED displays the charge status for the battery:

Color	Mode
Green flashing	Charging
Green	Battery fully charged

There are two receptacles in the front of the base station, used to hold light guides when they are not inserted in the handpiece.

5 - MODES USED

The **MINI L.E.D.** has three different modes. To select the mode you want to use, press the navigation button. The green LED indicator next to each symbol shows the mode selected.

"FAST CURING" mode: when the handpiece is removed from the base station, the device will be set on fast

curing program. This menu activates maximum light power (depending on the light guide) for 10 seconds.

"PULSE CURING" mode: Selected by means of the navigation key. This menu gives full power in a pulsed mode, which consists of an emission of 10 successive one second flashes with a rest period of 250ms between the flashes.

"STEPPED CURING" mode: Selected by means of the navigation key. "Stepped curing" (curing similar to the step menu in a halogen device), gives you:

- a soft start of 10 seconds
- full power for 10 seconds.

6 - EVERYDAY USE

- 1) The **MINI L.E.D.** is normally positioned on the base station. The first time you screw it onto the battery it will be activated in the fast curing mode. For further details, please refer to chapter 5.
- 2) Your **MINI L.E.D.** is ready to use once the menu has been selected.
- 3) Place the light guide as close as possible to the surface of the VLC material without touching it, as this could greatly damage the light guide and decrease its efficiency (staining)
- 4) Press the On/Off button briefly to start the curing cycle). The setting will be confirmed by an audible signal (beep).
- 5) At the end of the cycle, you will hear another beep. However, the cycle can be interrupted at any time simply by pressing the On/Off button.
- 6) If the device is not used, after three minutes it goes into Standby mode and all the indicator lights go out. Press any button to reactivate the device.

7 - CARE

Warning: For everyday maintenance of the **MINI L.E.D.**, make sure that the unit is not on its rest. Disconnect the rest from the **MINI L.E.D.** before using a disinfectant.

- 1) The **MINI L.E.D.** and its accessories can be cleaned and disinfected with ready-to-use cleaning and disinfecting wipes impregnated with alcohol, an

amphoteric and a biguamide (refer to manufacturer's instructions) for at least two minutes. Leave the product to take effect for at least 15 minutes. Use wipes that have CE marking or that are in compliance with any standard required by the national regulations.

Wipe dry using a disposable clean non-woven cloth until there are no traces of liquid.

Packaging: Keep sterilizable accessories in single-use sterilization sleeves or bags in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

Sterilize the light guide and the flexible optical shield individually in an autoclave before each patient, according to the following parameters:

- autoclave, compliant with standard EN 13060, class B
- sterilization temperature: 134 °C
- duration of the sterilization hold time: 18 minutes
- Pressure: 2 bars

Storage: Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is not conformable, recondition and re-sterilize in accordance with the defined protocol.

Disposition of the product: Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

- 2) **Before cleaning the handpiece, insert the protective cap in place of the light guide to prevent any liquid entering the handpiece.**
- 3) After each use, check thoroughly that no composite residue has adhered to the light guide. Remove any residue immediately and check that the surface of the light guide has not been altered. If damage is detected, replace the light guide, as the power of the curing light could be reduced significantly.
- 4) Under normal conditions of use, the power of the curing light does not vary if the battery is properly charged. Consequently there is no need to check the power as in the case of ordinary curing lights. However, if in doubt, check the power using the light guide.
- 5) The battery unscrews from the handpiece for replacement.
- 6) Remove the light guide to check that the LED is clean

and undamaged. Clean it with a dry air jet if necessary.

- 7) Never use an ultrasound bath for cleaning the **MINI L.E.D.** or its accessories.
- 4) The battery is located at the end of the handpiece. It can be removed by unscrewing.
- 5) The LED (located at the end of the handpiece) against which the light guide is fitted should regularly be checked for cleanness. If necessary, use a jet of perfectly dry air to clean it.
- 6) The light guide is not chemiclave safe. The conditions and sterilization cycles are according to the autoclave instructions. Steam sterilization: 134°C at 2 bars (200kPa) for 18min.
- 7) Disposal of the equipment and its accessories should be carried out by a collection centre for electronic equipment or by returning them to SATELEC's After Sales Service. In no case should the device and its accessories be disposed of by the user.

8 - TROUBLESHOOTING

If a problem appears, before calling your dealer's after sales service:

- 1) Check that the base station is connected properly to the mains to ensure that the battery recharges normally. If the LED is not lit up in spite of proper connection, contact the after sales service.
- 2) Check that the battery indicator shows that it is fully charged before pressing the On/Off button. If the LED is red, an automatic safety mechanism will stop the device from being activated. In this case, recharge your battery by replacing the handpiece on the base station, or just replace with a second fully charged battery.
- 3) Intensive use of the lamp will generally cause a rise in its external temperature that can be felt by the hand holding it. This is absolutely normal. If the temperature inside the **MINI L.E.D.** should rise too high, an automatic safety mechanism will stop the lamp, and the status indicator will flash red. Leave the device for a few minutes to cool down and then carry on as normal.
- 4) If the wattmeter control LED is red, check that the LED and the light guide are clean. Use a jet of dry air to remove any dust. If the problem persists or if the

light guide is damaged, the MINI L.E.D. must be returned to the after-sales service department.

9 - SPECIFICATIONS

Your dealer's after sales service is at your disposal to help with any technical problem.

Handpiece

Model: MINI L.E.D.
Weight: 165 g
Size: Ø23 x 200mm
Classification: Ordinary
Type B
Continuous service
IPX0

Power supply adapter

Input voltage: 110 - 240V AC
Frequency: 47 - 63 Hz
Output voltage: 12 V
Output current: 1,25 A
Classification: II
IP 40

Base station

Input voltage: 12 VDC
Protection: Fuse 2 A
Classification: Continuous service
IPX 0

Battery

Type: Lithium-Ion
Size: 90 x Ø21 mm
Capacity: 2000 mAh

Optical specifications

Wavelength: 420-480 nm
Intensity: 1250 to 2000 mW/cm² ± 10%
(depending on light guide)
class 2M laser

Temperatures

Operating: +10°C to +40°C
Storage: -20°C to +70°C

Humidity

Operating: 30% to 75%
Storage: 10% to 100% including condensation

10 - ELECTROMAGNETIC COMPATIBILITY

Warning: The charger power cord must be kept apart from those of any nearby devices.

MINI L.E.D. requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter 3.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the **MINI L.E.D.**. The recommended separation distances in this paragraph must therefore be complied with.

The **MINI L.E.D.** must not be used near or on top of another device.

If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the **MINI L.E.D.**.

10.1 - Electromagnetic emissions

The **MINI L.E.D.** is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the **MINI L.E.D.** is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1 Class B	MINI L.E.D. uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment. MINI L.E.D. uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.

10.2 - Magnetic and electromagnetic immunity

The **MINI L.E.D. ORTHO** is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the **MINI L.E.D. ORTHO** is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Magnetic field at 50Hz. IEC61000-4-8	3A/m	3A/m	The intensity of the magnetic field should be equivalent to that of a typical commercial or hospital environment (hospital, clinic).
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6KV contact ± 8KV air	± 6KV contact ± 8KV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2KV for power supply lines	± 2KV for power supply lines	Power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5.	± 1KV differential mode ± 2KV common mode	± 1KV differential mode ± 2KV common mode	Power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	<5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	Power quality should be that of a typical commercial or hospital environment. If the use of the MINI L.E.D. requires continued operation during a power cut, it is recommended that the product be powered from a separate power supply (UPS, etc.).

10.3 - Electromagnetic immunity / mobile radiofrequency equipment

The **MINI L.E.D.** is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	aCompliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the MINI L.E.D. (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz to 800MHz. $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below: 			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

10.4 - Recommended separation distances

The **MINI L.E.D.** is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled. The **MINI L.E.D.** user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the **MINI L.E.D.**, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in meters (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

10.5 - Cable lengths

Cables and accessories	Maximum length	Complies with:
charger power cord	<3m	RF emission, CISPR 1 - Class B Immunity to magnetic fields: IEC61000-4-8. Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3

11 - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 15), for information on how to proceed.

12 - LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.);
- repairs have been performed by persons not authorized by the manufacturer;
- the device has been used in an electrical installation which does not comply with current standards;
- the device has been used in a way which is not stipulated in this Manual;
- accessories other than those supplied by SATELEC have been used;
- the instructions in this document have not been followed.

The manufacturer reserves the right to modify the unit and/or the Operating Manual without notice.



13 - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

14 - SYMBOLS AND ABBREVIATIONS

SYMBOL	DEFINITION
	Alternating current
	Direct current
	Warning, please refer to the accompanying documentation
	Device emitting Class 2M laser rays
	Type B
	Class II

Note:

Technical personnel of the Satelec authorized dealer network can obtain from **ACTEON Group** on request all the information they need for repair of the parts of the curing light that ACTEON has identified as repairable.

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