# AIRSONIC ABSORBO BOX

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## GB Instruction of Use

### 1. Discription

The Airsonic Absorbo-Box is a universal transportable compact dust collector. Due to its small dimensions (36 x 38 x 22 cm) the Airsonic Absorbo-Box requires very little space in the dental laboratory. It only needs electric supply. Its integrated, powerful exhauster draws the abrasive dust into the filter element, and lets the clean air emit from the unit's backside. A bright fluorescent tube provides an optimal illumination of the work place. Elastic rubber sleeves on both sides of the Airsonic Absorbo-Box allow for reach into the working space. The Airsonic Absorbo-Box is switched on and off by a switch located on the left side of the unit.

### 2. Cleaning

It is essential to disconnect the unit from power supply before it is cleaned! The Airsonic Absorbo-Box should be cleaned at regular intervals. Heavy dust deposits in the filter reduce the exhauster's suction capacity. Clean the unit by opening the tilting lever on the right side of the unit's casing, remove the dust filter, and empty the Airsonic Absorbo-Box over a dustbin. Knock the filter carefully out against the inside of the dustbin. The dust filter should not be washed. The filter should be replaced in case it is perforated or ripped. The dust filter should also be replaced in case the airflow is still reduced after the filter is cleaned. The average economic life of the dust filter lies between 1 and 2 years.

### 3. Exchanging the Fluorescent Tube

It is essential to disconnect the unit from power supply prior to exchanging the fluorescent tube! Exchange the lamp by carefully pulling the fluorescent tube rightwards from the socket. Do not apply heavy pressure on the glass tube, as the tube may burst and injure your fingers. Use a compact OSRAM DULUX S/E 11W/840 tube for replacement. The average economic life of the tube is approximately 8,000 hours.

### 4. Exchanging the Fuse

It is essential to disconnect the unit from power supply prior to exchanging the fuse! In case the ventilator does not work even though the unit is turned on and connected to the power supply, either one or both fuses are defective. Once the power supply cable is disconnected the from the plug receptacle on the unit's left side, the fuse compartment can be pulled from the plug receptacle and the fuse can be removed. Use only 230 V, M 0.315 A micro fuses for replacement. Should the fuses fail again, the Airsonic Absorbo-Box is defective and should be sent to the manufacturer for repair.

### 5. Exchanging the Acrylic Window

It is essential to disconnect the unit from power supply prior to the exchange! In case the acrylic window is scratched through abrasive medium after extended use, it is easily exchanged by untightening the three fixation screws.

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