 HENRY SCHEIN®

IMPROVED FORMULA

EUROSIL MAX²

9008532 / 9008533 EUROSIL MAX² PUTTY

9008534 / 9008535 EUROSIL MAX² PUTTY, SOFT

9008536 EUROSIL MAX² CORRECTION

9008537 EUROSIL MAX² CATALYST LIQUID

9008538 EUROSIL MAX² CATALYST PASTE

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MANUEL D'UTILISATION
GEBRAUCHSANWEISUNG
ISTRUZIONI D'USO

MADE IN ITALY
Fabricado en Italia/Fabriqué en Italie
Hergestellt in Italien/Prodotto in Italia

Distributed by: (US only)
 HENRY SCHEIN INC.
135 DURYEA ROAD
Melville, NY 11747 USA

 HENRY SCHEIN U.K. HOLDINGS LTD.
MEDCARE HOUSE, CENTURION CLOSE
GILLINGHAM BUSINESS PARK
GILLINGHAM, ME8 0SB U.K.
www.henryscheinbrand.com

Rev. 2011/12

EUROSIL MAX²

C-SILICONES	PAGE 1 – 4	GB
SILICONAS POR CONDENSACIÓN	PÁGINA 5 – 8	ES
SILICONES PAR CONDENSATION	PAGE 9 - 12	FR
KONDENSATIONSVERNETZENDES SILIKON	SEITE 13 - 16	DE
SILICONI PER CONDENSAZIONE	PAGINA 17 – 20	IT
CONDENSATIESILICONE	PAGINA 21 – 24	NL
POLYSILOXANY	STRANA 25 – 28	CZ
YOĞ UNLAŞTIRILACAK SILICON	SAYFA 29 – 32	TR







Table 1

1	Eurosil Max² Putty Type 0 Putty Consistency	Eurosil Max² Putty Soft Type 0 Putty Soft Consistency	Eurosil Max² Correction Type 3 Light-bodied Consistency
2	30"	30"	30"
3	1'15	1'15	1'30
4	3'15"	3'15"	3'30"
5	4'30"	4'30"	5'00"
6	2% - 5%	2% - 5%	2% - 8%
7	> 98%	> 98%	> 99%
8	< -0,2%	< -0,2%	< -0,7%

GB

HYDROCOMPATIBLE CONDENSATION SILICONES (POLYSILOXANES) FOR HIGH-PRECISION IMPRESSIONS USE

Eurosil Max² Putty: Polysiloxane for the first impression characterised by initial flow in the working phase and high rigidity after hardening. Recommended for the two stage technique (double impression).

Eurosil Max² Putty Soft: Polysiloxane for the first impression characterised by initial flow in the working phase and adequate elasticity on removal from the oral cavity. Recommended for the two stage technique (double impression).

Eurosil Max² Light : Low viscosity, extremely flowable, hydro-compatible, polysiloxane. Recommended for the two stage technique (double impression).

IMPRESSION TRAY

All rigid types; Henry Schein stainless steel impression trays are particularly recommended.

EUROSIL MAX² Putty/ EUROSIL MAX² Putty SOFT - INSTRUCTIONS FOR USE

Dosage: Use the measuring spoon to dose the required quantity of Eurosil Max² and spread it out over your hand (note: for the correct dosage the putty should be leveled with the top of the measuring scoop). Press the rim of the measuring scoop onto the putty once for each measure dosed. For each measure of material, apply two strips of Eurosil Max² Catalyst Paste as long as the measuring scoop (4 cm.) When using Eurosil Max² Catalyst liquid, dose 5 – 6 drops of hardening liquid for each scoop of putty used.

A circular logo with the letters 'GB' inside, positioned at the top center of the page.
GB

Mixing: Fold and knead the putty energetically with your fingertips for about 30 sec., until the colour is even, without any stripes. Spread out the putty over the impression tray. The working time including mixing time is about 1 min. 15 sec.. Put into the mouth and wait until it is set, which will take about 3 min. 15 sec..

EUROSIL MAX² LIGHT - INSTRUCTIONS FOR USE

Dosage: Extrude the required quantity, observing the dosage scale on the mixing block. Then dose a quantity of Eurosil Max² Catalyst Gel equal to the length of the fluid. When using Eurosil Max² Catalyst Liquid, dose one drop of hardening liquid for each notch of light material.

Mixing: Mix energetically with the spatula, and then carefully spread the mixture out, pressing it over the mixing block to eliminate any air which may have been incorporated into the mixture. Repeat this operation until the mix is even in colour. Mixing time should ideally be about 30 sec. Apply the mixture to the impression tray with a spatula or a syringe for elastomers. Working time including mixing time is about 1 min. 30 sec. Put into the oral cavity and wait for the mix to harden, which will take about 3 min. 30 sec.

GB

TECHNICAL DATA

Table 1 gives the characteristics of the materials:

- 1 - ISO 4823 / ADA Spec. 19
- 2 - Mixing time
- 3 - Total working time*
- 4 - Time in mouth
- 5 - Setting Time*
- 6 - Strain in compression (min - max)
- 7 - Recovery from deformation
- 8 - Linear dimensional change 24 hours

* The times for clinical use (listed in **Table 1**) are intended from the start of mixing at a temperature of 23°C – 73°F. Mixing at higher temperatures reduces these times. Mixing at lower temperatures increases them.

WARNING / CAUTION

- Before making the corrective impression, remove all remaining traces of oral fluids by rinsing thoroughly and drying.
- Properly close the catalyst tube after use.
- Prevent catalyst from coming into contact with skin and eyes; wash well with soap and water in the event of accidental skin contact. In the event of eye contact, rinse thoroughly with running water and seek medical attention. If substance is ingested, seek medical attention immediately.
- Please avoid contact with garments.
- In sensitive subjects, polysiloxanes may cause irritation or other allergic reactions.

GB

CLEANING AND DISINFECTING THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected. Always follow the manufacturer's instructions.

CASTING MODELS

The impression must be cast 30 min. to 72 hours after hardening. The material is compatible with all quality commercially available chucks. We recommend Plaster type 4 die stone and type 3

GALVANISATION

The impressions may be galvanised with copper or silver.

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