HENRY SCHEIN® Z-DUPE

Addition-Vulcanizing Duplication Silicone.

Product description:

Addition-vulcanizing vinyl-polysiloxane in 3 different final hardnesses for universal use in the duplication technique. Z-DUPE distinguishes by a superior reproduction of details, constancy of shape and long shelf life. Due to its platinum hardener system the material is virtually shrinkless and it is not necessary to pour out the duplicating form immediately.

Dosage:

Z-DUPE is mixed in a volume ratio of 1:1 of the components catalyst (A) + base (B).

Processing:

Mixing of the two components can be effected either in a beaker or with the help of the numerous dosing and mixing devices available.

Processing time (incl. mixing):	6 ± 1 min.
Removal:	20 ± 3 min.
Final hardness:	15 ± 2 shore A 20 ± 2 shore A 28 ± 2 shore A

The mixing and processing time refer to a room temperature of 23 °C and a relative air humidity of 50 %.

Preparatory steps:

Almost every material can be duplicated with Z-DUPE, e. g. plaster models metals, etc. The master model should not be soaked with water, but should be dry and free of all insulating agents. To save material the base part

of the model can be embedded in the foamed plastic ring supplied with the duplication flask in order to cover only the teeth and jaw parts to be duplicated.

Before pouring out the duplicating form is to be sprayed with a surface-tension releasing agent!

This reduces the risk of air bubbles which could occur during pouring out. All kinds of dental plaster, model resins and embedding masses are suitable materials for the duplication model.

Further processing:

The best results regarding precision and homogeneity are achieved by means of the so called pressure duplicating technique. Put the just filled flask into a corresponding pressure pot. The master model may be deflasked directly after termination of the vulcanisation. The form is maintained and can be poured out several times.

If the duplicating form was produced under pressure it is decisive that the model is to be produced under pressure, too.

General information:

Z-DUPE shows no negative effects in physiological and toxicological respect. Chemical and thermal disinfection are possible.

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