Distributed by:

Distribuido por: Vertrieb durch:

Distribué par: Distribuito da:

**HENRY SCHEIN INC.**

Melville, NY 11747 USA

EC REP

Henry Schein U.K. Holdings Ltd.  
Gillingham ME8 0SB U.K.

 HENRY SCHEIN®

# NATURAL ELEGANCE SINGLE BOND

*UNIVERSAL LIGHT-CURING*

*DENTAL ADHESIVE*

## **INSTRUCTIONS FOR USE**

Agente de adhesion fotopolimerizable, multifunctional

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---

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Henry Schein U.K. Holdings Ltd.  
Gillingham ME8 0SB U.K.

# **GB - NATURAL ELEGANCE SINGLE COAT BOND**

**Light cured, self-etching adhesive system**

## **INSTRUCTIONS FOR USE**

**Henry Schein Natural Elegance Single Bond is a light-cured, multipurpose, one component adhesive agent for adhesive restoration techniques in dentistry. Use etching gels for etching of dentine and enamel before the application of Natural Elegance Single Bond is required.**

**Composition, Natural Elegance Single Bond contains:**

- Methacrylates • Polyalkenoate methacrylized • Amorphous silica

**Indications, Bonding adhesive restoration techniques**

Direct filling technique

- Adhesive bonding of composite materials to natural enamel and dentine
- Adhesive bonding of compomer materials
- Dentine sealing
- Sealing of sensitive tooth necks

Indirect restorative technique

- Adhesive bonding of ceramic and composite patterns with indicated cement to natural enamel and dentine.

Adhesive bonding of other dental materials

- Composite materials to pretreated ceramics • Composite material to composite material • Composite material to pretreated metals and amalgam

**Contra-indications**

If allergy exists to any of the components of Natural Elegance Single Bond. If the site cannot be isolated after enamel etching and during application and curing of Natural Elegance Single Bond. If oral hygiene is poor.

**Side effects**

Components of Natural Elegance Single Bond may lead to sensitization in patients with a special predisposition.

**Interactions with other agents**

Phenolic and other substances (e.g. zinc oxide-eugenol) which inhibit polymerization must not come into contact with Natural Elegance Single Bond. Contamination of etched enamel (e.g. by saliva) neutralizes the adhesive effect of Natural Elegance Single Bond.

**Application: A. Direct Filling Technique**

1. Cavity preparation

Prepare cavity according to the principles of the adhesive restoration technique. If indicated, apply hardsetting calcium hydroxide cement in areas near the pulp, using the spot technique. Use the rest of the dentinal surface for bonding.

## 2. Etching

Use etching gel according to the manufacturer's instruction. This is applied directly from the syringe onto all enamel and dentine surfaces to be used for bonding, and distributed evenly using a disposable brush. No traces of etchant should be visible on the enamel or on the dentine. Remove excess water from the surface with a cotton pellet or with a short spray from the air gun. Do not dry the surface completely, as the exposed collagen layer might collapse and reduce adhesion.

## 3. Application of Natural Elegance Single Bond

Apply Natural Elegance Single Bond directly from the syringe onto a disposable brush and massage into the cavity for 20 s. This time should not be shortened as it is required for adhesive diffusion into the collagen network. The massaging movement on the surface supports this process. Lightly air dry. Cure with halogen light for 30 s.

## 4. Composite-, Compomer filling

Place restorative according to the corresponding instructions for use.

### **B. Adhesive bonding of inlays/onlays, veneers and crowns**

#### 1. Preparation of the restoration

Preparation of the restoration is dependent on the material used and should be carried out according to the manufacturer's instructions.

#### 2. Cavity preparation

Remove temporary restoration, clean cavity, insert restoration to check the fit, then remove it.

#### 3. Etching

Use Etching according to manufacturer's instruction.

#### 4. Application of bonding agent

Apply Natural Elegance Single Bond on a brush, massage dentine and enamel for 20 s, lightly air dry, light cure for 30 s.

#### 5. Cementation

Cover all surfaces to be bonded evenly with a thin layer of selected cement, insert it into the cavity and firmly press into position. For optimum results wait 3 to 4 min, until the chemical setting process is finished. Cure each side for 60 s.

### **C. Adhesive bonding of other dental materials. C1 Fractured ceramic patterns, porcelain-fused-to-metal restorations and amalgam**

#### 1. Preparation of the restoration

Roughen metal surfaces and fractured ceramic surfaces using diamond abrasive tools (40  $\mu$ m). Etch ceramic surfaces with 5 % hydrofluoric acid (for extraoral use only) and apply a silane coupling agent. Etch surfaces (enamel, composite) to be repaired with 35 % phosphoric acid (etch gel for 30 s) (Etch Gel Henry Schein), rinse and dry. Optional: sandblast the metal. Clean surfaces with water or acetone, dry.

#### 2. Application of bonding agent

Apply Natural Elegance Single Bond directly to the metal or ceramic surface with a disposable brush. Lightly air dry, then light cure for 30 s.

### 3. Use of opaquer

Use an opaquer according to manufacturer's instruction after bond application.

### 4. Application of composite

Place restorative according to the instructions for use.

## **C2 Fractured composite restorations**

### 1. Preparation of the restoration

Roughen composite surfaces with a diamond abrasive tool (40 microns). Etch surfaces (enamel, composite) to be repaired with 35 % phosphoric acid (etch gel for 30 s) (Etch Gel Henry Schein), rinse and dry.

### 2. Bond application

Apply Natural Elegance Single Bond directly from the syringe onto a disposable brush and massage into the composite and/or tooth surface for 20 s. Cure with halogen light or LED light (like HS-LED or HS-CU light) for 30 s.

### 3. Application of composite

Place restorative according to the instructions for use.

## **D. Sealing of sensitive tooth necks**

1. Cleaning of the tooth: Clean the tooth with an oil-free polishing paste.

2. Conditioning: Proceed as described in A.2.

3. Bond application: Proceed as described in A.3.

We recommend a control-visit after 6 months.

## **Emergency measures**

Danger from swallowing can be avoided due to the small quantity applied (by means of a small brush). In case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In case of contact with the eyes, rinse thoroughly with water and seek medical advice.

## **Toxicological Effects**

The toxicity of the primer mixture that directly touches the vital dentine and – as intended – penetrates into the dental tubules can be compared to that of polycarboxylate or glassionomer-cement.

## **Note**

Only supplied to dentists and dental technicians. Keep out of the reach of children!

## **Shelf life and marking**

Expiration date and number are marked on the containers. Do not use after shelf life expiration date.

## **Storage**

Store at 4–23 °C / 39–73 °F. Avoid exposure to direct sunlight or other heat sources.

### **Caution**

Federal law restricts this device to sale by or on the order of a dentist.

### **Infection control**

Compliance with the following instructions ensures optimal hygienic conditions when using Natural Elegance Single Bond

- Use a new disposable brush for each patient. • Needles are for single use only.
- Close the syringe after use by replacing the cap and clean with a disinfectant.

**Date of issue:** June 2008

**For MSDS see [www.henryscheinbrand.com](http://www.henryscheinbrand.com)**

## **E - NATURAL ELEGANCE SINGLE COAT BOND**

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