 HENRY SCHEIN®

LAB VIBRATOR

VIBRADOR DE LABORATORIO

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HS-RÜTTLER

VIBRATORE DA LABORATORIO

LAB VIBRATOR

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MANUEL D'UTILISATION
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING

Distributed by:
Distribuido por: Vertrieb durch:
Distribué par: Distribuito da:

HENRY SCHEIN INC.
Melville, NY 11747 USA
Henry Schein U.K. Holdings Ltd.
Gillingham ME8 0SB U.K.

Users manual for LAB VIBRATOR

Vibrator with electronic power regulation

Table of Contents

GENERAL NOTES	1
DESCRIPTION OF THE EQUIPMENT AND ITS USE	1
TECHNICAL SPECIFICATIONS	2
IMPORTANT WARNINGS	2
INSTRUCTIONS FOR USE	2
MAINTENANCE	2
WARNINGS	3
MALFUNCTION CAUSES-REMEDIES	3
ELECTRICAL SCHEMES	3
SPARE PARTS	3

GENERAL NOTES :

Please read this booklet carefully. It contains important information about safe installation, use and maintenance. Keep this booklet in a convenient place for further reference.

- (1) Check the general condition of the unit after removing the packaging. When in doubt do not operate the unit and contact a qualified service repair person authorised by **HENRY SCHEIN**.
- (2) Before connection, check to see that the identification plate indicates the correct voltage for the electrical supply system. Installation must be carried out in compliance with all applicable safety standards and in accordance with the manufacturer's instructions. The manufacturer will not be held responsible for damage to objects and/or harm to people caused by improper installation.
The unit is safe only when connected correctly to ground via an appropriate device in accordance with safety standards.
This fundamental requirement must be satisfied. If in doubt, request a thorough inspection of the system by a professionally qualified technician.
The manufacturer will not be held responsible for damage caused by an improperly grounded system.
- (3) This unit should only be used as intended by the manufacturer. Any other use must be considered inappropriate and consequently dangerous. The manufacturer will not be held responsible for damage caused by improper, erroneous or unreasonable use.
- (4) Switch-off and unplug the unit from the electrical outlet before servicing.
- (5) Carefully follow cleaning instructions as in the manufacturer's manual.

DESCRIPTION OF THE EQUIPMENT AND ITS USE :

The vibrators of the **LAB VIBRATOR** series have been developed for vibrating and compacting investments, plasters or similar materials. This procedure improves the homogeneity of the product and eliminates air bubbles that could be generated while mixing and melting the product.

The **LAB VIBRATORS** are equipped with an easy to remove rubber plate to facilitate cleaning and a power regulator which allows continuous change to the vibration intensity from minimum to maximum. The **LAB VIBRATOR medium** model is also supplied with a two-way switch for selecting a lower and a higher vibration.

The **LAB VIBRATOR** series consists of 3 different models with different dimensions and powers:

LAB VIBRATOR small:	plate 16 x 15,5 cm;	max 100W
LAB VIBRATOR medium:	plate 22,5 x 14 cm;	max 200W
LAB VIBRATOR large:	plate 34 x 24 cm;	max 350W

TECHNICAL SPECIFICATIONS :

	LAB VIBRATOR S	LAB VIBRATOR M	LAB VIBRATOR L
Casing	Steel housing with baked-enamel finish	Steel housing with baked-enamel finish	Steel housing with baked-enamel finish
Dimensions (width x depth x height)	17 x 17 x 8,5 cm	24 x 15,5 x 10 cm	35 x 25 x 13 cm
Plate Dimensions	16 x 15,5 cm	22,5 x 14 cm	34 x 24 cm
Voltage (VERIFY THE INDICATION ON THE EQUIPMENT LABEL)	230Volt/50Hz or 120Volt/60Hz	230Volt/50Hz or 120Volt/60Hz	230Volt/50Hz or 120Volt/60Hz
Max absorbed power	100 Watt	200 Watt	350 Watt
Weight (approx.)	kg 3,3	kg 5,1	kg 8,7

Power regulator for the modification of the vibration intensity two way switch for lower and higher vibration (**on Lab Vibrator M only**)

IMPORTANT WARNINGS :

- (1) Check to see that the main voltage corresponds to the power requirements for the unit as specified on the identification label.
- (2) During the vibration the material on the vibrator could overflow. Avoid any spillage of material from the rubber plate in to the machinery as this could result in damage.

INSTRUCTIONS FOR USE :

- (1) Check to see that the main voltage corresponds to the power requirements for the unit as specified on the identification label. Plug cable into machine and power socket.
- (2) On Lab Vibrator medium select the desired vibration by using switch (17) while the machine is switched off.
- (3) To start vibrator use on/off switch (16).
- (4) Adjust the vibration level by means of the knob (15).
- (5) Place the ring with the material to be vibrated on the rubber place for the desired time.

MAINTENANCE :

- (1) The vibrators require minimal maintenance apart from occasional cleaning of the rubber plate. This should be carried out with the rubber plate removed from the vibrator.

WARNINGS :



- (1) If replacement of an internal part is required, disconnect the equipment from the power before removing the upper plate.

MALFUNCTION CAUSES-REMEDIES :**(1) The vibrator doesn't start:**

- 1.1 Check that the plug is correctly inserted in the sockets and that the electrical cord is not damaged.
 1.2 Check that the fuse is not burned out, otherwise replace it with one with the same characteristics.
 1.3 Let check from a to Have a qualified technician check that electricity is present at the main socket.

(2) The Vibrator lights but does not vibrate:

- 2.1 The power regulator (13) must be replaced.
 2.2 The coil (6) must be replaced.

ELECTRICAL SCHEMATICS :**LAB VIBRATOR small/large****LAB VIBRATOR medium****SPARE PARTS :**

Pos.	DESCRIPTION	Cod. LAB VIBRATOR S	Cod. LAB VIBRATOR M	Cod. LAB VIBRATOR L
3	Rubber plate	CGV010	CGV008	CGV009
4	Foot	CPP006	CPP004	CPP004
6	Coil	EBA003	EBA004	EBA005
7	Vibro stop	CEV001	CEV002	CEV002
11	Plug with fuse holder	EFS005	EFS005	EFS005
12	Electrical cord	EWA011	EWA011	EWA011
13	Electronic power regulator	EVC001	EVC001	EVC001
14	Min-Max panel	M10	M10	M10
15	Knob	EVM001	EVM001	EVM001
16	Switch	EID005	EID005	EID001
17	Two way switch (only for LAB VIBRATOR M)		EID004	
18	On-off indicator (only for LAB VIBRATOR L)			ELA003

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