# Henry Schein Vest™

Investment material for crown and bridgework for all alloys. Suitable for pressable ceramic as well as for use in speed and conventional technique

## Mixing ratio

Powder : Liquid 100 g : 24 - 26 ml  
Powder : Liquid 160 g : 39 - 41 ml

Concentrate mixtures for metal rings and ringless casting:	100g	160g
* 50 % Onlays, inlays, conicals crowns	12 liq + 12H <sub>2</sub> O	20 liq + 19H <sub>2</sub> O
* 65 % Gold castings, crowns and bridges	16 liq + 8H <sub>2</sub> O	25 liq + 14H <sub>2</sub> O
* 75 % Secondary parts and baked ceramics	18 liq + 6H <sub>2</sub> O	29 liq + 10H <sub>2</sub> O
* 100 % Non-ferrous metal alloys (non precious alloys)	24 liq	39 liq

The concentrate admixture figures are for guidance only - they depend on the type of alloy being used.

## Processing

- Measure Liquid/Powder and mix thoroughly for 30 sec
- Agitate under vacuum for 1 min
- Maintain the vacuum for another 30 sec

Adjust to the lowest agitation setting while pouring. After finishing, no further agitation is required.

## Processing time approx. 5 min

If using metal rings: line with ring liner as follows:

Muffle sizes 1 and 3: 1 insert  
Muffle size 6: 2 inserts

Casting without rings: use a packing which is smooth on the inside. Overlapping edges should be greased with vaseline.

## Speed casting alloys

### Preheating for 25 – 30 min.

After casting, the muffle should be place inside the oven at a temperature of 850-900°C.

## Important

When using ovens with floor heating, make sure there is sufficient clearance (approx. 1 cm) between the mould and the floor plate

## Preheating times

Muffle size	1	35 min
	3	45 min
	6	60 min

In each case at final temperature. For each additional ring increase the preheating time for another 10 minutes.

## Over night heating alloys

During the heating over night the liquid concentration should be reduced by 5% to 10%, according to alloy and application, as the full setting expansion is reached.

## Burnout

After hardening place the mould into the cold oven. At 290°C and 580°C the temperature should be maintained for 30 - 45 minutes, depending on the size and the number of the moulds. Hold the final temperature, recommended for the alloys, 30 - 45 minutes.

## Heating rate

app. 3 - 5 °C/min

## Maximum temperature: 1200 °C

Cast without delay.

Also suitable to pressable ceramics.

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Henry Schein U.K. Holdings Ltd.