HAEMOSTATIC AGENTS

• Instructions for use

FERRIC
SULFATE 20%

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INSTRUCTIONS FOR USE

Description:

Haemostatic Gel 20% Ferric Sulfate contains inert binding agents in a viscous, aqueous vehicle. The pH is ~1.0.

A. Haemostasis

1. A 1.2ml unit dose syringe is attached to the end of the Haemostatic bulk syringe and the agent is dispensed into small syringe. Caution: To prevent cross-contamination, a new unit dose syringe must be loaded for additional volumes.
2. Twist an infuser tip securely onto the 1.2ml syringe. Check evenness of flow before dispensing intra-orally.
3. Dispense solution while rubbing until adequate haemostasis is obtained. Firmly rub until all bleeding stops.
4. When haemostasis is obtained, clean preparation thoroughly with a firm air/water spray. If bleeding occurs with spray, repeat step 3 and give another firm air/water spray.

Precautions:

1. Haemostatics and or blood Haemostatic mixtures can prevent adequate adhesion and lead to staining microleakage under "direct bonded" restorations. Preparations must be thoroughly washed using a firm air/water spray and scrubbing with pumice. If any mucin, coagulum, or residual haemostatics are present they may contaminate the bonding interface and/or prevent polymerisation/set of resins.
2. Temporary cements and protein mucins can contaminate preparation surfaces for indirect "bonded/luted" restorations. Mineral astringents such as ferric sulfate can cause these materials to become more firmly attached to the tooth. Scouring the surface with pumice to ensure the preparation site is absolutely clean is imperative.
3. For definitive bonding and luting of restorations, thorough scouring and cleaning of the preparation is crucial to avoid micro-leakage and decreased bonding strength.
4. Haemostatic solutions are designed for intraoral use, however, check for evenness of flow before dispensing in the mouth.
5. Ferric sulfate haemostatics should not be locally mixed with epinephrine or any epinephrine-impregnated retraction cord since a temporary blue colouration may occur.
6. No practical benefit is gained by using ferric sulfates in conjunction with other agents; mixing ferric sulfate with aluminum chloride renders a solution more acidic than either one alone.
7. Please rinse Haemostatic thoroughly from application site using air/water spray. The set of polyether material may be affected if this is not done.
8. The presence of iron may cause a temporary discoloration of the gingival tissues. The patient should be reassured that this colour will disappear within 24-48 hours. The Haemostatic solutions are yellow-coloured but being water-soluble, they can be washed off and will not permanently discolour or stain enamel, dentin or gingiva. Some base mate-

rials stain on the surface but with no lasting effects. Some fabrics may stain with these haemostatics.

9. Do not use on patients with known allergy to ferric sulfate or on patients with chemical hypersensitivities.
10. Do not use self-etching primers or self-etching primer/bonding systems immediately after using haemostatics. If self-etching bonding agents are to be used after haemostasis, always scour surfaces with pumice and thoroughly wash before application. This is not necessary with conventional bonding systems that have a stand-alone etching step. It is also not necessary when using conventional glass ionomer, zinc phosphate, or similar cements following attainment of haemostasis.
11. Thoroughly wash Haemostatic from tissues before applying other chemistry
12. Use protective glasses when handling/loading haemostatic.

General Precautions:

1. For professional use only.
2. Review instructions, precautions, and MSDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight.
4. Avoid skin exposure to resins.
5. Isolate strong chemicals to area of treatment.
6. Confirm that patient has no known allergies to treatment materials.
7. Test flow of materials from syringe and tip before using intraorally.
8. Never force syringe plungers.
9. Constant high volume suction is necessary to prevent product from being swallowed.

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