 HENRY SCHEIN®

900-8181

DIGITAL BLOOD PRESSURE UNIT

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 0473

GENERAL

This instruction manual is intended to assist the user for safe and efficient operation of the automatic digital blood pressure monitor (hereinafter: device) model LD7. The device must be used in accordance with the procedures described in the manual. It is important to read and understand the entire manual, especially the section "Tips on taking blood pressure measurement". This device is intended for the non-invasive measurement of systolic and diastolic arterial blood pressure and pulse rate in adults (age 15 and above). Consult the physician if measurement is taken in children or persons with arrhythmia as errors may occur.

PRINCIPLE OF OPERATION

This device adopts the oscillometric technology with Fuzzy Algorithm measuring the arterial blood pressure and pulse rate. The cuff is wrapped around the arm and automatically inflated by the air pump. The device sensor catches weak fluctuation of the cuff pressure in the cuff produced by extension and contraction of the artery of the arm in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters from the mercury column, and is displayed in digital value. Annotation: This device cannot provide reasonable accuracy if used or stored in the temperature or humidity beyond the range stated in the section "SPECIFICATIONS" of this manual.


NEW TECHNOLOGIES USED

FUZZY algorithm takes into account the uniqueness of individual heartbeats, which provides higher accuracy of measurement.

TIPS ON TAKING BLOOD PRESSURE MEASUREMENT

1. Arterial blood pressure is subjected to sharp fluctuations. The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter.


Arterial blood pressure changes with atmospheric pressure and is affected considerably by many factors, e.g. physical loads, emotional state, stress, meals, etc. Medicines, drinking, smoking has a big effect on the individual level of blood pressure. When blood pressure is measured in hospital, the value is always higher than that at home. Tension can have a serious effect on certain groups of patients and is known medically as 'White coat effect'. Blood pressure will raise in low temperature, so it is better to take blood pressure measurement in room temperature (approximately 20°C). If device is stored in low temperature, bring up to room temperature for at least 1 hour before use or measurement could be inaccurate. Blood pressure does vary with age and individual, and it is recommended to write down the readings in blood pressure record daily and check with your doctor to find out what 'normal blood pressure' is for you.

2. Take measurement under doctor's instruction for patients with cardio-vascular diseases. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor!
3. Accurate measurement of blood pressure may be difficult in serious arteriosclerosis, weak pulse, or in patients with obvious fluctuation of heart contraction rhythm. Please consult qualified physician to interpret your blood pressure readings.
4. It is necessary to remain still and silent during measurement to get accurate readings. Measurement should be conducted in quiet environment at room temperature. Don't eat or smoke before a measurement. This device is supplied with the cuff which will fit arm sizes 22-42cm. Children and adults whose cuff size falls outside the range 22-42 cm should select another cuff size. Please contact the dealer to get the other cuff sizes.

ATTENTION:

Do not use cuffs other than the original cuff contained in this kit!

Arterial blood pressure classification by World Health Organization

5. Repeated measurements at 3 minute intervals are recommended to calculate the average for a more accurate measurement. Atherosclerosis patients should use longer intervals (10-15 minutes) as elasticity of patients' vessels decrease significantly in these diseases. 10-15 minutes interval is also applicable for patients suffering from diabetes for a long time. Per ESH (European Society of

Hypertension) recommendations 'At least two measurements at 1 min intervals are should be taken carefully at each visit', we provide the average of the latest three readings.

BATTERY INSTALLATION

1. Open the battery cover and install four 'AA' type batteries into the battery compartment as indicated. Make sure that the polarity is correct.
2. Close the battery compartment cover.


Replace batteries if W appears in the display or if no power when ON/OFF button is pressed.


Replace all batteries simultaneously, and don't use rechargeable batteries. If the device is to be unused for a long time, please take out the batteries. Do not leave worn batteries in the device. Life span of enclosed batteries may be short and should only be used to check that the device is operating correctly.

USE THE DEVICE WITH AC POWER Adaptor

In addition to the batteries, you can use AC power adaptor as the power supply. AC power adapter is optional for the device for sale. The connector of AC power adapter is located in right side of the device. Use only AC power adapter with below specified technical features Output voltage: $6V \pm 5\%$

- Max. output current: At least 600 mA
- Output plug polarity: <-> inner
- External diameter 5.5 ± 0.1 mm
- Internal diameter: 2.1 ± 0.1 mm
- Length: 10 ± 0.3 mm

CORRECT POSTURE FOR MEASUREMENT

1. Sit at a table and let the table support your arm as you take the measurement. Make sure that the cuff on the upper arm is at approximately the same level as the heart, and that the forearm is extended naturally on the table.
2. You may lie on your back and take measurement. Look at the ceiling, keep calm, and don't move your neck or body during the measurement. Make sure that the cuff on the upper arm is at approximately the same level as the heart.

FITTING THE CUFF

1. Insert the edge of the cuff approximately 5 centimeters into the D-ring as shown in figure.
2. Put the cuff on the left upper arm with the tube pointing to the direction of palm. If measurement on your left arm is difficult, you can use right arm for measurement. In this case, it is necessary to know that the readings may differ about 5-10 mmHg between left arm and right arm.
3. Wrap cuff around your upper arm with the lower edge of the cuff approximately 2-3 centimeters above the elbow. The mark "artery" must be over the artery of the arm.

4. Press the cuff to make sure that it is attached securely. The cuff should not be too tight or too loose. It should be easy to insert two fingers between the cuff and the upper arm.
5. The mark "INDEX" on the cuff must point to area "NORMAL" (22-32 cm). This means the cuff size is correct. If mark "INDEX" points to the area beyond area "NORMAL", please consult your dealer for another cuff size.
6. If it is difficult to secure the cuff evenly to the user's upper arm, a cone shape is also acceptable.
7. Clothes should not restrict the area of the upper arm, if sleeve is tight then remove clothing to allow accurate measurement.

SETTING THE DATE AND TIME


The function provides accurate measuring time for each measurement. Preset the date and time correctly before use as follows:

1. When the device is connected to power supply at first time, the display will show as Fig. 1.
2. Press and hold button "PRESET"
3. Press button "M1" or "M2" to add or subtract the number, and press button "START" for confirmation.
4. When the year setup is finished, the month number will flash automatically. Please follow the same instruction as above to set month, date and time.
5. Press button "START" to finish setup.


CARRY OUT A MEASUREMENT

1. Insert the tube plug into the air connector.
Before the measurement, take 3~5 times deep breaths and relax.
Don't talk or move your arm.
2. Press button "START " and all symbols will appear on display in 2 seconds as Fig. 3.


3. Two beeps will sound and '0 mmHg' will appear on the screen. Pump begins to inflate with display showing the reading of pressure. Generally the pressure will reach 190 mmHg as Fig. 4.

4. Pump stops inflating and pressure begins to decrease gradually, during which the user's blood pressure and pulse will be calculated as Fig. 5.

5. There will be a long beep once measurement is complete. The air in the cuff will quickly deflate and the blood pressure reading - pulse reading will show. The measuring time will also display in two screens and a reminder to record the reading will show " " flashing.

6. Press button "M1" or button "M2" to record the reading in corresponding memory. For example, if button "M2", is pressed, the display will show as Fig. 7. If the user does not press button, the reading won't be recorded.

7. If power supply is not switched off and the device is not used, it will return to standby mode after 3 minutes.

AUTOMATIC INFLATION

There are 4 given levels of inflation pressure for this device: 190mmHg, 230mmHg, 270mmHg and 300mmHg. When 190mmHg is not enough or movement of arm occurs, the device will automatically inflate to reasonable pressure level to ensure a successful measurement. It is not a fault.

RAPID DEFLATION DURING MEASUREMENT

If you feel unwell during measurement or want to stop for any reason. Press the On/Off button. The device will quickly release the air in cuff and the device will be returned to standby mode.

FUNCTION OF MEMORY

MEMORY RECALL


1. LD7 can store 60 sets of readings each in "M1" and "M2", and will automatically calculate the average value of the latest 3 readings for "M1" and "M2" respectively. When the memory is full (60 sets of readings are stored), the oldest reading will be replaced by a newer one. Memory remains even if power is removed.

2. Once measurement is finished or if in standby mode, the user can press button "M1" or button "M2" to recall memory.

Press button "M1" or "M2", the display will show the average value of the latest 3 readings as Fig. 8.


Press again, the display will show '01', which means the latest reading, then turns to another screen to show readings and measuring time as Fig. 9.

Press again, the display will show '02', which means the second to the latest reading...

MEMORY CLEARANCE

After a measurement is finished or when the device stands by, hold down button "M1" or "M2" for at least 5 seconds, the display will show "CLR", which means the stored reading for "M1" or "M2" is removed as Fig. 10.

ERROR AND LOW BATTERY INFORMATION

INDICATION	POSSIBLE REASON	CORRECTION METHODS
C Err	<ul style="list-style-type: none"> • The cuff is put on incorrectly or the tube plug is inserted too loosely • Movement of arm/hand or talking during measurement. • The cuff is not inflated to necessary pressure • Arrhythmia 	<ul style="list-style-type: none"> • Make sure that cuff is put on correctly and the tube plug is inserted tightly and repeat the measurement • Repeat the measurement the following completely recommendations of manual • Repeat the measurement with pumping cuff on 30-40 mmHg above expected systolic pressure • Consult your personal physician
–	The batteries are weak	Replace all 4 batteries with new ones

CARE, STORING, REPAIR AND RECYCLING

1. Protect this device from high moisture, direct sunlight, shock, solvent, alcohol and gasoline.
2. Remove batteries if device is stored for a long time, and keep the batteries away from the children.
3. Keep cuff away from sharp objects and do not extend or twist the cuff.
4. Clean device with soft or dry cloths only.
5. The cuff must be handled with care and the cover can be cleaned with a damp cloth. **WARNING:** Do not wash the inner bladder.
6. The device should be checked yearly by specialists to ensure correct operation. Please consult your dealer for more information.
7. Follow the local directive (WEE regulations) for disposal of batteries and device.

TROUBLESHOOTING

SYMPTOM	CHECKPOINT	REMEDY
No display when the ON/OFF button is pressed	<ul style="list-style-type: none"> • The batteries have run down • The polarity of battery is wrong • The contact of battery compartment is polluted 	<ul style="list-style-type: none"> • The batteries have run down • The polarity of battery is wrong • The contact of battery compartment is polluted
Inflation stops and re-inflate later	<ul style="list-style-type: none"> • The automatic inflation for ensuring correct measurement • Did you talk or move your arm (or hand) during measurement? 	<ul style="list-style-type: none"> • See "AUTOMATIC INFLATION" • Keep still and silent during the measurement

SYMPTOM	CHECKPOINT	REMEDY
The reading is extremely low or high.	<ul style="list-style-type: none"> • Is the cuff at the same level as the heart • Is the cuff wrapped correctly? • Did you strain your arm during measurement? 	<p>Make sure that your posture is right</p> <ul style="list-style-type: none"> ■ Wrap the cuff correctly • Relax during measurement • Keep still and silent during the measurement
Pulse rate is too low or too high	<ul style="list-style-type: none"> • Did you talk or move your arm (or hand) during measurement? • Did you talk or move your arm (or hand) during measurement? Did you make measurement right after exercise? 	<ul style="list-style-type: none"> • Keep still and silent during the measurement <p>Take measurement again after resting for more than 5 minutes</p>
The batteries are run down soon	<ul style="list-style-type: none"> • Faulty batteries are used 	<p>Suggest to use alkaline batteries of known manufacturers</p>
The device is automatically turned off	<p>It is the result of automatically turn off system</p>	<ul style="list-style-type: none"> • This is to save the power consumption of the device, and it is not a fault.

WARRANTY OBLIGATIONS

1. The Monitor has a 36 months warranty from date of purchase excluding the cuff. The cuff has a 12 month warranty.
2. The warranty obligations are prescribed by warranty certificate for buyer.
3. The addresses of organizations for guarantee maintenance are present in the warranty certificate.

SPECIFICATIONS

Model	LD7
Size	150 (U x 130 (W)x 50(H) mm
Weight	Approximately 350 g without batteries
Measuring method	Oscillometry
Measuring range	40 to 260 mmHg (blood pressure) 40 to 160 beats/minute (pulse rate)
Measuring accuracy	± 3 mmHg for systolic and diastolic pressure $\pm 5\%$ of the reading for the pulse rate
Inflation	Automatic by the pump
Rapid deflation	Automatic electronic valve
Batteries	4"AA"x1,5V
Adapter	Optional component, 6V, 600 mA

LD7 DIGITAL BLOOD PRESSURE MONITOR GUARANTEE AND REPAIR SERVICE

Your digital blood pressure monitor is covered by the best service and guarantee policy in the industry. The digital blood pressure monitor is guaranteed against manufacturing faults and defective materials for a period of 1 year. If any defective materials or manufacturing faults are found during the term of the guarantee, Henry Schein will make the necessary repairs free of charge on receiving the faulty device. The guarantee does not cover accidental damage or evident misuse of the device. In the event of any claim, please contact our distributor.

Date of purchase

Distributor

important notice! THIS IS A COMPLEX TECHNICAL DEVICE. READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

EN1060-1:1996/EN1060-3:1996

ANSI/AAMI SP10

IEC6060-1

IEC6060-1-2

The manufacturer reserves the right to make technical changes without notice in the interest of progress.

Prior notices will not be given in case of any amendments within this manual. The mentioned trademarks and names are owned by the corresponding companies.

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