Distributed by (in US only)

HENRY SCHEIN INC.
135 DURYEA ROAD
Melville, NY 11747 USA

EC REP

HENRY SCHEIN U.K. HOLDINGS LTD.
MEDCARE HOUSE, CENTURION CLOSE
GILLINGHAM BUSINESS PARK
GILLINGHAM, ME8 0SB U.K.
www.henryscheinbrand.com

Rev.2013/08

Maxima® by  HENRY SCHEIN®

NATURAL ELEGANCE 7th Generation

Single Bond

Universal adhesive • self-etching, single component light-cured

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

Maxima[®]by **NATURAL ELEGANCE 7TH GENERATION**

IS A SINGLE COMPONENT LIGHT-CURED, SELF-ETCHING ADHESIVE FOR ADHESIVE RESTORATION TECHNIQUES.

INSTRUCTIONS FOR USE

Composition

Natural Elegance 7th Generation contains:

Methacrylates

Photo initiators

Ethanol

Water

Indications for use

Adhesive bonding for use with adhesive restorative techniques:

- Adhesive bond of composite (e.g. Maxima Nanocomposite from Henry Schein) and compomer materials (e.g. Compomere from Henry Schein) to enamel and dentin
- Adhesive bond of ceramic and composite restorations (e.g. Maxima Nanocomposite from Henry Schein) to enamel and dentin
- Adhesive bond of composite material to ceramics, composite material (e.g. Maxima Nanocomposite from Henry Schein), metal and amalgam (e.g. Dental Amalgam from Henry Schein)
- Dentin sealing

Contraindications

Hypersensitivity due to any of the elements within Natural Elegance 7th Generation. Inadequate oral hygiene. If the working area cannot be kept completely dry during application.

Side effects

Elements of Natural Elegance 7th Generation may lead to sensitivity or an allergic reaction for patients with a predisposed condition.

Interactions

Phenolics and other substances (e.g. zinc oxide eugenol) may not come in contact with Natural Elegance 7th Generation, as they will inhibit polymerisation.

Application

Note

To prevent the bond from curing in the dispensing well, Natural Elegance 7th Generation should be immediately covered with the black cover after each application.

1. Conventional filling technique

1.1 Cavity preparation

Cavity preparations should be prepared according to the principles of the adhesive restoration technique. If indicated, apply hard-setting calcium hydroxide cement (e.g. Henry Schein Calcium Hydroxide) in areas near the pulp, using the spot over-capping technique. Prepare the enamel using a finishing bur. In cases with enamel that is not prepared, etch it using phosphoric acid (e.g. Henry Schein Etch Gel) or cleanse it using pumice (e.g. Henry Schein Pumice). Thoroughly rinse and dry the enamel afterwards. For hygienic reasons, it is recommended to use a dental dam (e.g. 900-2651 HS Latex Dental Dam or 900-3055 HS Non-Latex Dental Dam).

1.2 Application

Shake the bottle well before use and dispense one drop into the dispensing well. Massage Natural Elegance 7th Generation using a disposable brush for 20 s onto the dentin and enamel. Gently dry for 5 s using oil-free compressed air.

Light-cure for 10 s (e.g. Henry Schein Halogen- or LED Polymerisation units with a light output of >800 mW/cm²). If using a curing lamp with a light output less than this, be sure to increase the length of the curing time accordingly.

1.3 Composite and compomer fillings

Apply restorative material according to the manufacturer`s instructions for use.

2. Adhesive bonding with inlays/onlays, veneers and crowns

The application of Natural Elegance 7th Generation is only possible under the condition that sufficient light access (according to the manufacturer`s instructions for use) can be achieved through the complete restoration.

2.1 Preparation of the restoration

The preparation of the restoration is dependent on the material used and according to the manufacturer`s instructions for use.

2.2 Cavity Preparation

Proceed as described in point 1.1.

2.3 Application

Proceed as described in point 1.2. Avoid over-application of the adhesive as this could impact the fit of the restoration.

2.4 Cementing

According to the manufacturer`s instructions for use.

3. Adhesive bonding with other surfaces

3.1 Preparation

3.1.1 Ceramics

Intraoral application:

Roughen ceramic surface using a diamond bur or sandblasting, etch with 35 % phosphoric acid, rinse, and dry using oil-free compressed air.

Extraoral application:

Etch ceramic surface using 5 % HF-acid, rinse, and dry using oil-free compressed air.

Optional:

Ceramic surfaces can also be silanised (according to the manufacturer's instructions for use).

3.1.2 Metal and amalgam

Roughen metal surfaces using a diamond bur or sandblasting. Clean using steam jet (extraoral) or acetone; and dry using oil-free compressed air.

3.1.3 Composite

Roughen composite surfaces using a diamond bur or sandblasting, etch with phosphoric acid, rinse and dry using oil-free compressed air.

3.2 Application

Shake the bottle well before use and dispense one drop into the dispensing well. Afterwards close the bottle. Massage Natural Elegance 7th Generation using a disposable brush for 20 s onto the surface. Gently dry for 5 s using oil-free compressed air.

Light-cure for 10 s (e.g. Henry Schein Halogen- or LED Polymerisation units with a light output of >800 mW/cm²). If using a curing lamp with a light output less than this, be sure to increase the length of the curing time accordingly.

3.3 Composite and compomer fillings

Apply restorative material according to the manufacturer's instructions for use.

4. Dentin sealing

Proceed as described in points 1.1 and 1.2.

Emergency measures

In case of direct contact with oral mucosa, it is sufficient to rinse area using tap water. In case of contact with eyes, rinse eyes thoroughly using tap water. Seek immediate professional medical attention from an ophthalmologist or eye care specialist.

Notes

For use only by dentists and dental technicians, or upon their request. Keep out of the reach of children!
Only recycle completely empty containers.

Shelf life and marking

Expiry date and number are marked on the immediate containers and packaging.

Do not use after the expiry date.

Storage at temperature 2–8 °C / 35,6–46 °F.

Keep out of direct sunlight or near sources with warm temperatures.

Infection control

Compliance with the following recommendations will help ensure optimal hygienic conditions, when using Natural Elegance 7th Generation:

- Only use disposable brush for one patient.
- After each use, close the bottle with the cap, and disinfect using conventional methods.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Date of issue

08-2013

Maxima[®]*by* **NATURAL ELEGANCE 7TH GENERATION**

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Composicion

Natural Elegance 7th Generation se compone de:

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