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Tandheelkundige thermoformende film

**Eigenschappen**

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1,0 / 1,5 / 2,0 / 3,0 / 4,0 / 5,0 mm

**Colours**

nature-transparent – 3 mm blanks in various

colours available.

**Values of the material**

Density:	0.95 g/cm <sup>3</sup>
Softening temperature:	approx. 43 °C
Tensile strength:	33 MPa
Tear strength:	> 18 N/mm <sup>2</sup>
Tear elongation:	850 %
Hardness:	82 Shore A
Water absorption:	< 0.2 %

**Application**

Duplication of models, soft splints, positioners, mouthguards or bleaching.

**Processing**

All Dreve thermoforming units are suitable for the processing of DrufoSoft® blanks (Drufoformat scan, DrufoSmart scan, DrufoSmart, Drufoformat SQ / TE,

Vacformat U / P, Vacformat 2000). For plastification the blank will be placed into the plate reception and with the help of the radiant heater brought to the necessary temperature. The times vary between 1:30 and 2:45 minutes according to the unit type and the surrounding temperature. After end of the heating time the material hangs down → complete the thermoforming process. Before thermoforming plaster models should be watered and insulated. Depending on the material thickness the cooling time under pressure / vacuum is 3–8:30 minutes. After the cooling time DrufoSoft® can be cut by means of scissors or a scalpel. DrufoSoft® is suitable for „additive thermoforming“. Through process of lamination various foils can be connected (**only with pre-forming unit under pressure**)

**Information on health**

DrufoSoft® is marked with a CE label and therefore fulfils the general requirement of the medical device guideline.

**Keep away from children!**

**Only for dental use!**

**Storage conditions**

At room temperature, i. e. 18–28 °C.

**Durability:** see date of expiry.

DrufoSoft® is not to be used when expiry date has passed.

**This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to**

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