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Addition-vulcanising vinyl-polysiloxane with a final hardness of 22 Shore A for all dental duplications. Easy and precise to mix in a ratio of 1:1.

Product description

Dublisil® 22 plus is an universal duplication material for all kinds of duplications. It distinguishes by excellent mechanical properties, as e. g. increased tear-resistance and tear-growth-resistance. Due to the practical way of mixing and the very thin-flowing initial consistency absolutely precise duplication models can be produced. Dublisil® 22 plus-forms are nearly shrinkless and have a high form consistency as well as a long shelf life due to the used platinum hardener system. Thanks to the improved mechanical properties an unproblematic demoulding is guaranteed for. Negative forms made of Dublisil® 22 plus can be filled easily and for several times without losing precision. The surface of the duplication model is smooth and hard without after-treatment.

Preparatory measures

Almost every material, e. g. plaster models, model resins and metals can be duplicated with Dublisil® 22 plus **Important:** The master model must not be watered, but should be dry and free of any insulating agents etc. For an economical consumption of duplication material we recommend to use the Dublisil® flask (REF D3805) with a variably adjustable interior as well as the foam space maintainer (REF D38076) for blocking out the base parts which are unimportant for the duplication.

Dosage

Dublisil® 22 plus is mixed in the ratio 1:1 of the components A (catalyst / white) and B (basis / red).

Processing

The two single components can be mixed manually in a mixing beaker or automatically by means of a dosage and mixing unit. The most rational and clean method is mixing in the Dosper evo (REF D5050). At the touch of a button Dublisil® 22 plus is dosed, mixed and dispensed fully automatically without air bubbles. If the material is mixed manually, the mass is to be filled into the flask in a preferably thin stream. In case of the Dosper evo this stream is preset by the unit. The distance

This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence

between unit and flask should be approx. 20 cm.

Further processing

The best results regarding precision and homogeneity are achieved by the so-called pressure duplication technique. For this the just filled flask is to be put into a pressure pot (e.g. Polymax). In case that the duplication form is produced under pressure it is important that also the model is produced under the same pressure. After the vulcanisation has been completed the master model can be deflasked directly. The period of time between the deflasking of the master model and the production of the duplicate should to achieve full resilience of the silicone should be at least 30 minutes. Before filling the duplication form can be moistened with a relaxant. This leads to a lower surface tension and improves the flow properties of the model material. All dental plasters, model resins and embedding materials are suitable.

General aspects

Dublisil® 22 plus is unproblematic in physiological and toxicological respects. A chemical and thermal disinfection is possible. Wear protective clothing! Silicone impurities on clothes cannot be cleaned – also not chemically. Allergic reactions to product components cannot be excluded in case of sensitive persons. Vulcanisation can be disturbed by latex gloves, hand creams, cleaning agents etc. The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.

General information on this product is based on detailed application practice and intensive research work. Due to further development we reserve the right to modify the product.

Keep away from children!

Storage conditions

Store at room temperature, i. e. 18–28 °C.

Batch number / use-by date

The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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