ACRYL ORTHO (GB)

Directions for use

Colour: clear

Sprinkling the powder

Carefully apply a separating agent to the model. A high-quality alginate separator should be used. In case that the model is very dry, briefly soak it in water prior to applying the separator. Sprinkle the powder onto the model with the plastic bottle. Moisten the powder with Acryl Ortho monomer. Once the entire liquid has been soaked up and an homogeneous dough has been formed, apply the next layer. Build-up the desired form layer by layer. Once the entire appliance has been formed, soak the surface once more. The upper layer will quickly lose its shine and the model can be placed in the pressure vessel and the acrylic polymerized.

Polymerization under pressure

The water in the pressure vessel should be at 37°C (100°F). The acrylic must be fully immersed in water. A pressure of 2 - 4 atmospheres must be applied and maintained for 20 minutes in order to polymerize the

acrylic. When Acryl Ortho auto polymerizing acrylic is processed in this manner its density is so high that its mechanical properties correspond to those of heat curing acrylics. Acryl Ortho does not tend to form bubbles or porosities and exhibits a high precision of fit after being cured in a pressure vessel.

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Henry Schein U.K. Holdings Ltd.