Instructions for Assembly

1. Insert slide in front opening of dispenser gun.
2. Simultaneously raise metal bar on back of dispenser gun and push slide through opening.
3. Raise top hinge of dispenser gun and place MaxiTemp™ cartridge into channel.
4. Secure cartridge by snapping plastic hinge down into place.
5. Immediately prior to use of the material, remove cartridge cap and dispose.
6. Attach cannula tip by aligning notched side to cartridge and turning clockwise until fastened.

NOTE: It is important not to remove mixing tip after dispensing the material.

The used tip functions as a self-seal of the cartridge until the next use.

Instructions for Use

1. Fabricate a pre-impresion of the quadrant before any preparation is done to the teeth.
2. Insert MaxiTemp™ into the dispenser gun as per assembly directions.
3. Affix a new mixing cannula when you are ready to begin a new application of the MaxiTemp™ paste.
NOTE: It is important not to remove mixing tip after dispensing the material. The function of the used tip is to provide a self-seal of the cartridge until the next use.
4. Depress your dispenser gun handle to express MaxiTemp™ into your pre-impresion.
 - a. Express a small amount of the material prior to extrusion into your impression to insure a proper mix.
 - b. Working time to extrude paste and seat the filled matrix onto the prepared teeth is 45 seconds.
 - c. Begin placement of MaxiTemp™ into your impression on the occlusal area to eliminate the formation of air bubbles.
5. Seat filled impression onto the prepared teeth and wipe away any overflow with a plastic tool.
6. After 60 seconds, detach temporary from the mouth and allow to "bench cure" for 5 minutes.
7. When curing is complete, your temporary may be trimmed and polished with low speed acrylic burs and/or high speed finishing diamond burs.

NOTE: Trimming the MaxiTemp™ will produce composite dust. The use of safety glasses is recommended.

8. The oxygen inhibition layer is easily removable with ethyl alcohol.
9. Use a non-eugenol temporary cement in placing the finished crown or bridge into the mouth.

Repairing Broken or Cracked Temporaries

1. Roughen the break area and add a few undercuts to proximate sections.
2. Cover break area with a new mix of MaxiTemp™ and press pieces together for 60 seconds.
3. Excess material is detached after total curing time of 6 to 7 minutes.

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