Maxima® by HENRY SCHEIN®

Glass Fibre Posts & Reamers RADIOPAQUE & LIGHT TRANSMITTING

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GB

Instructions for use:

The Henry Schein fibre post is a highly radiopaque (350% Al.) and light conducting ($\geq 5\text{mW}$) fibre post. Thanks to its cylindro-conical structure, it is made to fit to the natural anatomy of the root canal.

Indications for use:

In case of insufficient residual tooth structure a fibre post is needed to support the coronal restoration.

Step-by-step:

- Select the correct post size according to the anatomical situation.
- Remove the root filling using a Gates Glidden Drill (rotation speed 800-1200 rpm).
- Calibrate the root canal using the corresponding Henry Schein Reamer (rotation speed 1000-2000 rpm.)
- Check post for fit, trim if necessary using a diamond disc and clean with alcohol.
- Apply etchant, primer, bonding agent according to manufacturer's instructions, or recommended self-etch bonding system, using a micro-applier.
- Apply the composite cement in the root canal and on the post, insert post into root canal and light cure.
- Commence core build-up immediately using a composite core build-up material.

Recommendation:

Henry Schein recommends disinfection of the rotary instruments using appropriate cleaning agent and sterilisation in an autoclave (134°) before each use.

Rubber dam is recommended whenever possible.

Safety and liability:

Henry Schein declines any liability for possible damage due to:

- Use of instruments not compatible with the system.
- Non-observation of the instructions for use.
- Usage outside of the field of application stated in the instructions for use.

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