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## Indications for use

Thermoforming blank for bleaching- and pharmaceutical splints.

## Product description

BioBleach® soft was especially designed as a splint material for home-bleaching. Also after the thermo-forming process the special formulation of the 1.0 mm strong EVA-material offers a consistent soft and stable splint. The usage for several bleaching processes is possible without any problems.

## Preparative measures

Before thermo-forming the blank a space maintainer as a reservoir for the bleaching material has to be applied on the labial areas of the teeth. For this exist two possibilities:

a) With a light-curing easy flowing but thixotropic resin material (Cenit L) the areas at the plaster model which are supposed to be bleached are consistently covered with a brush and afterwards get cured.

or

b) The whole model is covered with the thermo-forming blank Drufofen W 0.8 mm. Following the part of the model to be bleached (incl. a slightly overlapping of the cutting edges) is cut out of the blank by means of a scissors and put back at the model again.

**IMPORTANT:** In order to avoid any contact between the musoca and the bleaching material the space maintainer should end approx. 1 mm in front of the marginal gingiva.

## Processing instruction

Insert the BioBleach® soft blank in the thermo-forming unit and heat it. The right plastification is achieved when the blanks sags for approx. 10–15 mm. For Dreve-thermoforming units please take the heating times from he following chart:

	Drufofen SQ/TE	Vacfofen 2000	Vacfofen U/P
<b>BioBleach® soft</b>	1:05 min.	0:50 min	1:05 min.

After having achieved the right plastification drag the blank directly over the model and form it with vacuum or pressure.

## Processing instruction

The detaching of the splint from the cooled down thermo-forming blank is effected by means of a scissors or a hot scalpel. Due to the good elasticity of the splint material the marginal limit of the blank can outreach the gingival margin. This extends the vent edge of the splint and therefore guarantees the safe position of the bleaching material in the reservoir. In order to avoid irritations of the gingiva and the cheek the cut edges of the thermoforming splint should carefully be rounded off with a rotating instrument (e. g. UltraTrimm light) and smoothed.

## Important advice

• The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.

• The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

## Keep away from children!

**Only for dental use!**

## Storage conditions

At room temperature, i. e. 18–28 °C.

**Durability:** see date of expiry.

BioBleach® soft is not to be used when expiry date has passed.

## Form of delivery

**D4300** **BioBleach® soft** 20 blanks

This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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