**Remarques importantes**

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Soft denture relining resin.

## Indication for use

To be used as long-term soft relining material for removable dentures by direct or indirect method.

## Product description

**MegaBase®** is an autopolymerising silicone and is suitable both as a direct relining material and as a laboratory material. Apart from the relining material, **MegaBase®**-Set includes an adhesive for making a safe bond between silicone and PMMA resins, as well as a self-curing, two-component lacquer for final gloss coating.

## Fields of application

- Permanent soft relining
- Relief of permanent pressure spots
- Improvement of the suction effect (flabby ridge)
- Protective cover following implants
- Cover for sharp-edged bone parts

## Technical data

### Silicone

Processing time: ≥ 2:00 min

### Direct method / Setting time

In the mouth: 4:00 min  
At room temperature: 20:00 min  
In the water bath 35–40 °C (95–104 °F): 10:00 min

### Indirect method / Setting time

At room temperature: 20:00 min  
In the water bath 35–40 °C (95–104 °F): 12:00 min  
Final hardness: 35 ± 2 Shore A

### Adhesive

Setting time: 1:30 min

### MegaBase® 2K Lack

Working time: ≥ 2:15 min  
Setting time: 20:00 min  
room temperature: 10:00 min  
in the water bath 35–40 °C (95–104 °F): 10:00 min

The processing time refers to a room temperature of 23°C ± 1°C (73,4°F ± 1,8°F) and a relative air humidity of 50 %. Lower temperatures prolong, higher temperatures accelerate the processing and setting times. Technical data according to DIN EN ISO 4823.

## Processing

### 1. Preparatory measures

The prosthesis needs to be cleaned of residues of saliva and then dried. The total basal area is to be relined, as well as approx. 2–3 mm of the vestibule / buccal areas must be ground down by approx. 1 mm. Old relinings must be completely removed and the edges of the prosthesis smoothed. For a clean transition from the silicone to the denture base resin, a circumferential, rectangular stage should be incorporated. Impurities can significantly weaken the silicone / plastic composite. Therefore, once the preparatory measures have been completed, it is recommended that you clean and degrease the entire prosthesis with an alcoholic cleaning agent.


### 2. Processing of the adhesive

The adhesive is applied uniformly thinly and comprehensively to all smoothed, clean areas of the prosthesis where **MegaBase®** is to adhere, twice at an interval of about 30 sec. The drying time is about 1:30 min. Only continue working once the adhesive is completely dry.

### 3. Preparation of the cartridge system:

Before attaching the mixing tip, dispense material until it flows regularly out of both openings. Then fasten the mixing tip lock.

**Important:** After use, the mixing tip has to remain on the double cartridge until next needed.

## 4. Relining

### a) Direct method

The prosthesis is prepared as described in points 1 and 2. Then cover the entire base coat with a minimum thickness of 2 mm, starting from the edges of the prosthesis.

The processing time of **MegaBase®** is about 2:00 min. Apply all the material, and immediately incorporate it in the oral cavity. Until the processing time is complete, the patient carries out the necessary functional movements at normal occlusion. After 4 min. in the

mouth, the curing of the **MegaBase®** has progressed so far that the prosthesis can be removed. Before reworking (see point 5) with scissors and UltraTrimm the soft remaining relining should be put in a water bath at a temperature of 40–45 °C for a further 10 min so that the optimum physical properties are achieved. Alternatively, the prosthesis may also be stored for about 20 min at room temperature.

### b) Indirect method

The prosthesis prepared by the dentist from impression compound is inserted into a commercially available rebasing device or in a cuvette, and then processed as described after removal as at points 1 and 2. To avoid partial bonding between silicone and plaster, the model can be sealed with an alginate (e. g. **Isolat film**). By means of a dispensing unit **MegaBase®** is then applied everywhere basally as well as on the prepared vestibule / oral areas. To avoid air pockets, **MegaBase®** should be applied approx. 1 mm more than is necessary for the relining. The coated prosthesis is now reset in the rebasing device, and it is slowly and evenly built up pressure. The curing of the **MegaBase®** under layer takes around 20 min at room temperature. In order to accelerate the curing process, the relining can also be put in a pressure cooker or oven at around 40 °C for about 12 min.

## 5. Treatment

Excess material can be quickly and easily cut off with a scalpel or scissors. We recommend the use of rotating grinding wheels with different grain sizes (e. g. **UltraTrimm medium** or silicone abrasive heads) for subsequent smoothing. Please make sure the handpiece is at low speed (max. 8000 rpm).

## 6. Lacquering

See the **MegaBase® 2K Lack** instructions.

## Disinfecting

Use a 2 % glutaraldehyde or another appropriate solution per manufacturer's instructions. Disinfection procedure will not alter the soft denture relining resin.

## Important advice:

- Due to technical filling reasons there might be an unequal transportation of the two components of new cartridges. For this reason it is recommended to dispose the initially mixed pea-sized amount.
- Never apply using a mixing tip with contents that are curing or have cured. This can cause leaks in the cartridge system.
- It cannot be ruled out that some sensitive individuals may experience allergic reactions to product components.
- Curing may be impaired by the use of latex gloves, hand creams and cleansing agents.
- Remove all residues of registration material from the patient's mouth.
- Mixed material is chemically persistent. Spots onto clothing cannot be removed.
- This product is intended solely for the application described above and may only be used by qualified and trained personnel.
- We reserve the right to make changes to the product for development purposes.

## Keep away from children!

## Storage conditions:

Store at room temperature, i.e. 18 °C to 28 °C.

**Durability:** See expiry date.

Do not use **MegaBase®** after the expiry date.

## Adhesive

**Contains:** ethyl acetate, methyl methacrylate

## Hazard statements

H225 Highly flammable liquid and vapour. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. H336 May cause drowsiness or dizziness.

## Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P241 Use explosion-proof electrical / ventilating / lighting equipment. P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P405 Store locked up. P501 Dispose of contents / container in accordance with local / regional / national / international regulations.Vorschriften.

## Forms of delivery:

**REF** D5038223

**MegaBase®** Starter Kit  
1 x 50 ml double cartridge  
**MegaBase®** silicone  
1 x 50 ml double cartridge  
**MegaBase® 2K Lack**  
2 x 6 mixing tips  
1 x 10 ml adhesive  
2 x UltraTrimm medium  
incl. mandrell  
1 x grinding cap no. 2

**REF** D5038223NR

**MegaBase®** silicone  
2 x 50 ml double cartridges

**REF** D38233

**MegaBase® 2K Lack**  
2 x 50 ml double cartridge  
12 x mixing tips

**REF** D5038220NA

**MegaBase®** Adhäsiv  
10 ml bottle

**This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from inferior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.**

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