Características

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Field of application

Drufolen H is especially suitable for the production of hard-flexible cast form parts for the c & b technique and for the template technique at c & b temporaries. Drufolen H does not bond with acrylate and burns without residues.

Processing

Regarding the heating times of Dreve thermo-forming units please refer to the scale below:

Drufolen H	Drufomat SQ/TE	Vacfomat U/P	Vacfomat 2000
0,60 mm	1:20	1:20	1:00
0,80 mm	1:30	1:50	1:10

Before the thermo-forming process the plaster models should be watered and insulated with an insulating agent on alginate basis (Isolat film). By means of vacuum or pressure the plastified plate is formed on the object. The cooling down phase under pressure / vacuum should at least correspond to the respective heating time. Drufolen H can be reworked by means of a scissors. For the smoothing of the incisions special flexible grinding wheels (UltraTrim) are suitable. The processing with a hand adapter is also possible.

For the production of templates for the c & b temporaries the blank is thermo-formed on the situation model and cut in the areas of the of the ungrounded adjacent teeth which are located in the area which should be cared.

After the preparation the blank template can be repositioned in the patient's mouth resp. on the plaster model. The exact fixation is effected above the remaining teeth. The hollow space which results in the area of the preparation is filled with a temporary c & b resin (e. g. Fixtemp) and put back into the mouth / onto the insulated model. Please refer to the instruction of use of the corresponding resin materials for any further information regarding the processing.

Information on health

Drufolen H is marked with a CE label and therefore fulfils the general requirement of the medical device guideline.

Lot number / Expiry date

The lot number and the expiry date are indicated on each Drufolen W packaging. In case of claims, please always indicate the lot number of the product. Do not use the product after the expiry date.

Keep away from children! Only for dental use!

Storage conditions

At room temperature, i. e. 18 °C–28 °C.

Durability: see date of expiry.

Drufolen H is not to be used when expiry date has passed.

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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