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Tandheelkundige thermoformende film

Eigenschappen

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Water absorption: < 0.2 %

Application

Duplication of models, soft splints, positioners, mouthguards or bleaching.

Processing

All Dreve thermoforming units are suitable for the processing of **DrufoSoft® pro** blanks (**DrufoMat scan**, **DrufoSmart scan**, **DrufoSmart**, **DrufoMat SQ / TE**, **VacfoMat U / P**,

VacfoMat 2000). For plastification the blank will be placed into the plate reception and with the help of the radiant heater brought to the necessary temperature. The times vary between 1:45 and 2:15 minutes according to the unit type and the surrounding temperature. After end of the heating time the material hangs down → complete the thermoforming process. Before thermoforming plaster models should be watered and insulated. The cooling time under pressure / vacuum is at least 7 minutes. After the cooling time **DrufoSoft® pro** can be cut by means of scissors or a scalpel. **DrufoSoft® pro** is suitable for „additive thermoforming“. Through process of lamination various foils can be connected (**only with pre-forming unit under pressure**)

Information on health

DrufoSoft® pro is marked with a CE label and therefore fulfils the general requirement of the medical device guideline.

Keep away from children!

Only for dental use!

Storage conditions

At room temperature, i. e. 18–28 °C.

Durability: see date of expiry.

DrufoSoft® pro is not to be used when expiry date has passed.

suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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