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
## Working Instructions (EN)

### Odontosil 2K Lack

**Product description:**  
**Odontosil 2K Lack** is a self-curing sealing lacquer based on addition-cross-linked polyvinylsiloxanes for coating dental positioners based on A-silicones. **Odontosil 2K Lack** forms a thin, highly transparent and scratch-resistant surface. It is also absolutely free of solvents and thus completely odourless.

<b>Technical data:</b>	
Processing time 23 °C (73 °F)	2:15 min
Curing time 23 °C (73 °F; including processing time)	15:00 min
Final polymerisation (storage in water at 35-40 °C (95-104 °F))	15:00 min
Final hardness	25 ± 2 Shore A


The values specified for processing and curing times are valid for an ambient temperature of 23°C and a relative humidity of 50%. Lower temperatures prolong and higher temperatures shorten these times.

**Important advice:**  
 For technical reasons, there might be an unequal distribution of the two components in a new cartridge. For this reason, we recommend disposing of the first pea-sized amount. Never try to squeeze material through a cannula containing material that is curing or has cured. This may lead to leakages in the cartridge system and consequently to the discharge of unmixed and uncured individual components. Check the cartridge for uniform distribution before every use.

## Processing:

1. Prepare the silicone and smoothen the surface (e. g. with **UltraTrimm medium / light**).
2. Degrease the surface and let it dry in air for at least 3 min.
3. **Odontosil 2K Lack** is provided in convenient double cartridges, the mixing tips included in the scope of delivery allow an optimum dosage and mixture of the 2-component material.
4. Apply the required quantity either directly onto the surface to be coated or fill the material into a clean vessel. Then spread it thinly and evenly over the silicone with a brush. Possible spillage can be easily removed with a brush. The processing time is 2:15 minutes at room temperature (23 °C; 73 °F).
- 5a. The curing will begin independently once the two components have been mixed. The curing time is 15 minutes at room temperature.
- 5b. After curing the covered surface has to be brought in a water quench with a temp. of 35-40 °C (95-104 °F) for final polymerisation and best possible bonding towards the base silicone for 15 min.
6. Processing and curing times are prolonged by temperatures below 23 °C and shortened by temperatures > 23 °C.

**Important advice:**  
 For technical reasons, there might be an unequal distribution of the two components in a new cartridge. For this reason, we recommend disposing of the first pea-sized amount. Never try to squeeze material through a cannula containing material that is curing or has cured. This may lead to leakages in the cartridge system and consequently to the discharge of unmixed and uncured individual components. Check the cartridge for uniform distribution before every use.

**Keep away from children! Only for use in the area of application specified above!**

**Storage conditions:**  
 Store at room temperature, i.e. 18 °C–28 °C. **Shelf life:** See expiry date. **Odontosil 2K Lack** is not to be used after its expiry date has passed.

**Form of delivery:**  
**Art.-Nr. D38234**

**Odontosil 2K Lack**  
 2 x 50 ml double cartridges inc. 12 mixing cannulas, blue

**Accessories:**  
**Item No. D5007**  
**D4970X3**  
**D49701X3**

This product has been developed for dental use and must be applied according to the instructions. The manufacturer shall assume no liability in the event of damage resulting from misuse or improper application. Furthermore, the user is responsible for checking that the product is suitable for the intended purpose, especially if this purpose is not mentioned in the instructions.

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