# Bluephase® G4

The intelligent curing light

**EN Instructions for Use**

Page 2

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Pagina 56

**Appendix**

Page 110

**DE Gebrauchsinformation**

Seite 20

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Página 74

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Page 38

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Página 92

C US

**CE 0123**

**Rx ONLY**

For dental use only!  
Made in Austria

**Manufacturer**  
Ivoclar Vivadent AG  
9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

In USA distributed by  
Ivoclar Vivadent Inc.,  
175 Pineview Drive,  
Amherst, NY 14228, USA

**ivoclar**  
**vivadent**<sup>®</sup>  
*clinical*

*Dear Customer*

*Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The curing light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased Bluephase® G4.*

*Bluephase® G4 is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.*

*These Instructions for Use will help you safely start up the device, make full use of its capabilities and ensure a long service life.*

*Should you have any further questions, please do not hesitate to contact us (see addresses on the reverse page).*


*Your Ivoclar Vivadent Team*

# Table of Contents

<b>1. Product Overview</b>	<b>4</b>
1.1 List of parts	
1.2 Indicators on the charging base	
1.3 Indicators on the handpiece	
1.4 Operating the light	
<b>2. Safety</b>	<b>7</b>
2.1 Intended use	
2.2 Indications	
2.3 Signs & symbols	
2.4 Safety notes	
2.5 Contraindications	
<b>3. Start-Up</b>	<b>10</b>
<b>4. Operation</b>	<b>13</b>
<b>5. Maintenance and Cleaning</b>	<b>16</b>
<b>6. What if ... ?</b>	<b>18</b>
<b>7. Warranty /Procedure in Case of Repair</b>	<b>19</b>
<b>8. Product Specifications</b>	<b>19</b>

# 1. Product Overview

## 1.1 List of parts

## 1.2 Indicators on the charging base

Charging base without integrated radiometer:


- Indicator is black = Battery is charged
- Indicator flashes blue: battery is charging

Charging base with integrated radiometer and Click & Cure function (optional):

- Indicator is black = Battery is charged
- Battery shown on the indicator: battery is charging
- Light intensity shown on the indicator: measurement is carried out

### 1.3 Indicators on the handpiece

## 1.4 Operating the light

### Activating/deactivating polyvision

The polyvision function is activated by long pressing (>2 s) the program or time selection button (see 4. Operation). To confirm that polyvision is activated, the handpiece will beep and vibrate for a short time and the polyvision indicator will light up. Renewed long pressing of the program or time selection button results in the polyvision function being deactivated; the handpiece will not vibrate.

With the handpiece switched on, the current charging status is shown on the handpiece as follows:

• **No indicators lit up on the handpiece:**

**Battery sufficiently charged**

Curing capacity of minimum 20 minutes in the High Power program.

• **Battery symbol on the handpiece is flashing orange:**

**Battery weak**

Time/intensity can still be set and a polymerization time of approximately 3 minutes in the High Power program is left. Place the light into the charging base as soon as possible!

• **Battery symbol on the handpiece is flashing orange and a red "x" is shown:**

**Battery completely discharged**

The light can no longer be called up and the curing time can no longer be set. However, the handpiece can be used in the Click & Cure corded operation (only if using a charging base with an integrated radiometer).

## 2. Safety

### 2.1 Intended use

Bluephase G4 is an LED curing light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately in the oral cavity of patients. The intended place of application is in the dental practice, medical practice or in the hospital by the dentist or dental assistant. The intended use also includes the observation of the notes and regulations in these Instructions for Use. The curing light must only be operated by trained dental personnel.

### 2.2 Indications

With its "Polywave®" LED with broadband spectrum, Bluephase G4 is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385 – 515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporaries, as well as luting materials for brackets and indirect restorations, such as ceramic inlays.

### 2.3 Signs and symbols

Contraindication

*Symbols on the curing light*

Double insulation  
(device complies with safety class II)

Protection against electrical shock  
(BF type apparatus)

Observe Instructions for Use

Observe Instructions for Use

Caution

The curing light must not be disposed of in the normal domestic waste. Information regarding disposal of the light can be found on the respective national Ivoclar Vivadent homepage.

Recyclable

AC voltage

DC voltage

### 2.4 Safety notes

Bluephase G4 is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives IEC 60601-1-2 (EN60601-1-2) Edition 3.0, as well as the 93/42/EEC Medical Device Directive. The curing light complies with the relevant EU regulations.

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

### 2.5 Contraindications

Materials, the polymerization of which is activated outside the wavelength range of 385 – 515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not charge or use the appliance near flammable or combustible substances.

Never use without light guide.

The use of a light guide other than the one provided in the delivery form is not admissible.

The use of this device close to other equipment or stacked with it should be avoided because the correct function can be disrupted. If such use is unavoidable, the devices need to be monitored and checked for correct function.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Caution – The use of controls or adjustment devices or performing procedures other than those specified herein may result in hazardous radiation exposure.

### Warning

This unit should not be used near flammable anaesthetics or mixtures of flammable anaesthetics with air, oxygen or nitric oxide.

### Usage and liability

- Bluephase G4 must only be employed for the intended use. Any other uses are contraindicated. Do not touch defective, open devices. Liability cannot be accepted for damage resulting from misuse or failure to observe the Instructions for Use.
- The user is responsible for testing Bluephase G4 for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of the curing light at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent. The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light guide is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the handpiece during operation.

### Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

If the battery or power pack are used separately, e.g. during start-up or Click & Cure corded operation, contact with patients or third parties must be prevented. Do not touch the exposed contacts of the connection plug (power pack).

### Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the device is visibly damaged or no longer works correctly. A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source. Ensure that the device can be quickly and easily disconnected at any time.

### Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. To optimize user safety, the device has been equipped with intelligent anti-glare protection. For this, the "polyvision" function must be activated (see 4. Operation). With the polyvision function activated, the Bluephase G4 automatically senses if the handpiece is outside of the mouth and automatically switches off the light if it has been activated inadvertently.

Therefore, using the supplied anti-glare cones is recommended. Individuals who are sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to its light and wear protective orange goggles that absorb light below a wavelength of 515 nm. The same is true for patients.

### Battery

Caution: Use only original spare parts, particularly Ivoclar Vivadent batteries and charging bases. Do not short circuit battery. Do not touch the battery contacts. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that lithium-polymer batteries may react with explosion, fire and smoke development if handled improperly or mechanically damaged. Damaged lithium-polymer batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire and smoke development are toxic and corrosive. Do not touch leaking batteries with bare hands. In case of accidental contact with the eyes or skin, immediately wash with copious amounts of water. Avoid inhalation of fumes. In case of indisposition, see a physician immediately. Remove electrolyte residue from surfaces by washing/wiping with a moist cloth. Wash contaminated pieces of clothing immediately.

### Heat development

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

Generally, the stipulated curing times, particularly in areas near the pulp (adhesives: 10 seconds), must be observed.

Uninterrupted curing times of more than 20 seconds on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. Polymerize indirect restorations at intermittent intervals of 20 seconds or use external cooling with an air stream. The instructions regarding curing programs and curing times must be observed (see Selecting the curing program and the curing time). In addition, the light emission window must be placed exactly on the material to be cured at all times (e.g. by holding it in place using a finger).

**After several curing cycles on the same tooth, there is a risk that the pulpa suffers damage caused by the increased temperature!**

### 3. Start-Up

Check the delivery for completeness and any possible transportation damage (see List of parts). If parts are damaged or missing, contact your Ivoclar Vivadent representative.

#### Charging base without integrated radiometer

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply. Connect the power cord with the power supply. Make sure that the power cord is easily accessible at any time and can be easily disconnected from the power supply.

#### Charging base with integrated radiometer

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply.

Slide the connection plug of the power pack into the socket on the underside of the charging base. Tilt it slightly and apply slight pressure until you hear and feel it snap into place. Place the charging base on a suitable, flat table top.

Connect the power cord with the power supply and the power pack. Make sure that the power cord is easily accessible at all times and can easily be disconnected from the power supply. The charging base will briefly display "Bluephase G4" as ticker text on the screen.

### Handpiece

Unpack the handpiece from the packaging and detach the light guide by pulling it out. Then clean the handpiece and light guide (see Maintenance and Cleaning). After cleaning, reinsert the light guide.

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient (see Maintenance and Cleaning). Make sure to fit the protective sleeve snugly to the light guide. Then, attach the anti-glare cone or anti-glare shield to the light guide.

### Battery

We recommend fully charging the battery before the first use! If the battery is fully charged, it features a curing capacity of approximately 20 minutes. Slide the battery straight into the handpiece until you hear and feel it click into place.

Gently place the handpiece in the corresponding rest in the charging base without using any force. If a hygiene sleeve is used, please remove it before you charge the battery. If possible, use the light always with a fully charged battery. This will prolong the service life. It is therefore recommended to place the handpiece into the charging base after each patient. If the battery is fully discharged, the charging time is 2 hours.

Since the battery is an expendable part, it has to be replaced after its typical life cycle has expired after approximately 2.5 years. See battery label for the age of the battery.


### Battery: charging status

The current charging status is displayed on the handpiece as described on page 6.

**Click & Cure corded operation  
(only in charging base with integrated radiometer)**

If a charging base with an integrated radiometer is used, Bluephase G4 can be used in corded operation at any time, but particularly when the battery is completely empty.

For corded operation, release the battery by pressing on the release button on the lower part of the handpiece and then pull the battery out of the handpiece completely.

Then remove the power pack from the underside of the charging base. Do not pull on the power cord.

Insert the connection plug straight into the handpiece until you hear and feel it click into place.

During corded operation, the charging base cannot charge the battery, since it is not connected to a power source.

A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power socket.

**Measuring the light intensity  
(only in charging base with integrated radiometer)**

The integrated radiometer allows the light intensity ( $\text{mW}/\text{cm}^2$ ) to be measured easily and quickly while the charging base is connected.

To measure the light intensity, place the tip of the light guide, without protective cover, flush into the marked recess on the upper side of the charging base. Then activate the light and read the value displayed on the screen. If the light intensity is below  $400 \text{ mW}/\text{cm}^2$ , the screen will display "LOW".

## 4. Operation

Disinfect contaminated surfaces of the curing light as well as light guides and anti-glare cones before each use. Additionally, the light guide can be sterilized using the autoclaves intended for this purpose (see chapter Maintenance and Cleaning). Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For that purpose, check the light guide for contamination and damage as well as the light intensity at regular intervals (see paragraph Measuring the light intensity).

### Selecting the curing program and time

Bluephase G4 is equipped with the following 3 selectable curing times and 2 curing programs for the different indications. Use the time/program selection button to adjust the desired curing time and hence the specified light intensity.

#### H\* (High Power Program), 1,200 mW/cm<sup>2</sup> ± 10%:


The following curing times can be selected in the High Power program: 10, 15 or 20 seconds.

#### PRE (PreCure Program), 950 mW/cm<sup>2</sup> ± 10%:

The PreCure program is used to tack cure light-curing, adhesive luting composites (e. g. Variolink Esthetic) in order to facilitate the removal of excess material. The curing time of the PreCure program is preset to 2 seconds and cannot be altered.

**The PreCure program must not be used for conventional light-curing!**

Observe the Instructions for Use of the material applied when selecting the curing time and intensity. The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of max. 2 mm. Generally, these recommendations apply to situations where the emission window of the light guide is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to

be extended accordingly. For instance, if the distance to the material is approx. 11 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

- 1) The information provided herein applies to the 10 mm light guide supplied in the delivery form.
- 2) The information regarding heat development and burn hazards must be taken into consideration (see Safety notes).

Curing times		HIGH POWER 1,200 mW/cm <sup>2</sup> ± 10%
Restorative materials	<b>Composites</b> • 2 mm <sup>1)</sup> IPS Empress Direct / IPS Empress Direct Flow Tetric EvoCeram / Tetric EvoFlow Tetric / Tetric Basic White	10 seconds
	• 4 mm <sup>2)</sup> Tetric EvoCeram Bulk Fill / Tetric EvoFlow Bulk Fill Tetric Basic White	10 seconds
	Heliomolar / Heliomolar HB / Heliomolar Flow All conventional composites <sup>1)</sup>	15 seconds
	<b>Compomer</b> <sup>3)</sup> Compoglass F / Compoglass Flow	20 seconds
Indirect restorations / cementation materials	Variolink Esthetic LC <sup>4)</sup> / Variolink Esthetic DC <sup>5)</sup>	per mm ceramic: 10 seconds per surface
	Multilink Automix <sup>5)</sup> / SpeedCEM Plus <sup>5)</sup>	per mm ceramic: 20 seconds per surface
Adhesives	Adhese Universal ExcITE F / ExcITE F DSC Heliobond / Syntac	10 seconds
Temporary materials	Telio CS Link Telio CS Inlay/Onlay Telio Stains Telio Add-On Flow	10 seconds per surface 10 seconds 10 seconds 15 seconds
Miscellaneous	Helioseal / Helioseal F / Helioseal Clear	20 seconds
	Heliosit Orthodontic	10 seconds
	IPS Empress Direct Color IPS Empress Direct Opaque MultiCore Flow / Multicore HB	20 seconds 20 seconds 20 seconds

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

3) Applies to a maximum layer thickness of 3 mm

4) Applies to light-curing

5) Applies to dual-curing

### **Polyvision - Automatic assistant with intelligent anti-glare protection**

Bluephase G4 features a fully automated "polyvision" assistant for safe polymerization. This feature enables the curing light to automatically sense if the handpiece is inadvertently moved from its position during the curing procedure. To prevent any associated reduction in the amount of energy being transferred, the device will start to vibrate to alert the user to the improper use and, if necessary, will automatically extend the curing time by 20%. If the handpiece is significantly shifted away from its initial position (e.g. out of the oral cavity), the light will automatically switch off so that the curing process can be restarted and carried out correctly.

In addition, polyvision also acts as intelligent anti-glare protection. Polyvision ensures that the curing light cannot be called on in open space. The light can only be activated once the light guide is positioned directly over the material to be polymerized. This protects the operator and patient from being blinded by the light. If a protective sleeve is used, this function is not applicable.

If you do not want to use the assistant, you can deactivate it at any time by long pressing (> 2 seconds) the time or program selection button. The symbol on the handpiece (see 1.3 Indicators on the handpiece) disappears.

### **Measuring the light intensity (optional, only in charging base with radiometer)**

We recommend checking the light intensity emitted by the curing light at regular intervals using the integrated radiometer in order to be sure that the materials are adequately cured and the composite fillings are of a high and lasting quality (see 3. Start-Up).

### **Cure Memory function**

The last settings used, together with the combination of curing program and curing time, are automatically saved.

### **Start/Stop**

The light is switched on by means of the start/stop button. It is recommended that the emission window of the light guide is placed directly on the material to be polymerized. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start/stop button again.

### **Acoustic signals**

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Curing time and program change
- Inserting battery
- Error message (if anti-glare protection is activated and the curing process is aborted)

### **Light intensity**

The light intensity is maintained at a consistent level during operation. If the supplied 10 mm light guide is used, the light intensity has been calibrated to 1,200 mW/cm<sup>2</sup> ± 10% in the High Power program. If another light guide than the one supplied is used, it directly influences the emitted light intensity.

In a light guide with parallel walls (10 mm), the diameter of the light entry and the light emission window is the same. When using focussing light guides (e.g. Pin-Point light guide 6>2 mm), the diameter of the light entry is larger than that of the light emission window. The incident blue light is thus bundled to a smaller area. In this way, the emitted light intensity is increased. Pin-Point light guides are suitable for spot-on polymerization, e.g. to fix veneers before the removal of excess. For complete polymerization, the light guide must be changed.

## 5. Maintenance and Cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light guide. Disinfect contaminated surfaces of the device and anti-glare cones (FD 366/Dürr Dental, Incidin Liquid/Ecolab). Sterilize the light guide before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

### Cleaning the housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soapy solution.

### Pre-treating the light guide

Before cleaning and/or disinfecting the light guide, pretreat it. This applies to both automated and manual cleaning and disinfection:

- Remove substantial contamination immediately after use or 2 hours after that at the latest. For this purpose, thoroughly rinse the light guide under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to prevent the adherence of blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula, if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

### Manually cleaning and disinfecting the light guide

For manual cleaning, immerse the light guide in cleaning solution. Make sure that the light guide is entirely submerged in the solution and allow it to soak for the recommended time (ultrasonic cleaning or careful brushing with a soft brush may enhance the effectiveness of the procedure). We recommend using a neutral enzymatic detergent.

When cleaning and disinfecting, please make sure that the agents used are free of:

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solution (the maximum admissible pH value is 8.5)
- oxidizing agent (e.g. hydrogen peroxide)

Afterwards, remove the light guide from the solution and thoroughly rinse it under running, germ-free water (for at least 10 seconds).

For disinfecting the light guide, immerse it in disinfectant solution; make sure that the light guide is entirely submerged in the solution. We recommend using a o-phthalaldehyde-containing disinfectant.

After the light guide has been disinfected, remove it from the solution and rinse it again thoroughly under running germ-free water (for at least 10 seconds). Then dry the light guide with a clean towel.

Please observe the instructions given by the manufacturer of the detergent and disinfectant.

### **Mechanical cleaning and disinfection of the light guide (washer-disinfector)**

As an alternative, cleaning and disinfection can be carried out mechanically. Information on the validated procedure can be obtained from Ivoclar Vivadent AG.

### **Sterilization of the light guide**

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. The sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Dry the sterilized light guide using either the special drying program of your steam autoclave or hot air. The light guide has been tested for up to 200 sterilization cycles.

### **Checking the light guide**

After that, check the light guide for damage. Hold it against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light guide with a new one. If you can still see signs of contamination on the light guide, the cleaning and disinfecting procedure must be repeated.

### **Disposal**

The curing light must not be disposed of as normal household waste. Dispose unserviceable batteries and curing lights according to the corresponding legal requirements in your country. Batteries must not be incinerated.

## 6. What if ...?

Indicator	Causes	Error rectification
<p data-bbox="91 220 210 258">Red "x" lights up</p> 	<p data-bbox="241 220 456 242">The device is overheated.</p> <hr/> <p data-bbox="241 338 483 386">Electronic component of the handpiece is defective.</p>	<p data-bbox="582 220 1028 290">Allow the device to cool down and try again after a certain time. If the error persists, please contact your dealer or your local Service Centre.</p> <hr/> <p data-bbox="582 338 1021 411">Remove and reinsert the battery. If the error persists, please contact your dealer or your local Service Centre.</p>
<p data-bbox="91 462 210 517">Red "x" and battery symbol light up</p> 	<p data-bbox="241 462 360 485">Battery empty</p> <hr/> <p data-bbox="241 561 423 584">Battery contacts dirty</p>	<p data-bbox="582 462 1016 485">Place the device in the charging base and charge it.</p> <hr/> <p data-bbox="582 561 983 584">Remove battery and clean the battery contacts.</p>
<p data-bbox="91 711 215 775">The charging base is not illuminated during charging</p>	<ul data-bbox="241 705 506 778" style="list-style-type: none"> <li>- Power pack not connected or defective</li> <li>- Battery fully charged</li> </ul>	<p data-bbox="582 705 1021 778">Check if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord.</p>

## 7. Warranty / Procedure in Case of Repair

The warranty period for Bluephase G4 is 3 years from the date of purchase (battery: 1 year). Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- Damage resulting from improper handling, especially incorrectly stored batteries (see Product Specifications: Transportation and storage conditions).
- Damage to components resulting from wear under standard operating conditions (e.g. battery).
- Damage resulting from external influences, e.g. blows, drop to the floor.
- Damage resulting from incorrect set-up or installation.
- Damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- Damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, battery, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Centre. In case of a defect that cannot be rectified, please contact your dealer or your local Service Centre (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

## 8. Product Specifications

Light source	Ivoclar Vivadent Polywave® LED
Wavelength range	385–515 nm
Light intensity	High Power program: 1,200 mW/cm <sup>2</sup> ± 10 % PreCure program: 950 mW/cm <sup>2</sup> ± 10 %
Operation	3 min on / 7 min off (intermittent)
Light guide	10 mm, autoclavable
Signal transmitter	acoustic every 10 seconds and every time the start/stop button or the time / program selection button is activated or the anti-glare protection is enabled or the curing process is aborted
Dimensions of the hand-piece (without light guide)	L = 170 mm, B = 30 mm, H = 30 mm
Weight of the handpiece	135 g (including battery and light guide)
Operating voltage handpiece	3.7 VDC with battery 5 VDC with power pack
Operating voltage of the charging base without radiometer	100–240 VAC, 50–60 Hz, max 0.1 A
Operating voltage of the charging base with integrated radiometer	5 VDC
Power supply	Input: 100–240 VAC, 50–60 Hz max 1 A Output: 5 VDC / 3 A Manufacturer: EDAC POWER ELEC. Type: EM1024B2
Operating conditions	Temperature +10 °C to +35 °C Relative humidity 30 % to 75 % Ambient pressure 700 hPa to 1060 hPa
Dimensions of the charging base	D = 110 mm, H = 55 mm
Weight of the charging base without radiometer	155 g
Weight of the charging base with radiometer	145 g
Charging time	Approx. 2 hours (with the battery empty)
Power supply of the handpiece	Li-Ionen battery (approx. 20 min. with a new, fully charged battery in the High Power program)
Transportation and storage conditions	Temperature –20 °C to +60 °C Relative humidity 10 % to 75 % Ambient pressure 500 hPa to 1060 hPa The curing light has to be stored in closed, roofed rooms and must not be exposed to severe jarring. Battery: – Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15–30 °C / 59–86 °F – Keep the battery charged and store no longer than 6 months.
Delivery form	1 Charging base with power cord and power pack 1 Handpiece 1 Handpiece support 1 Light guide 10 mm 1 Anti-glare shield 3 Anti-glare cones 1 Pack of sleeves (1x 50 pcs) 1 Instructions for Use

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# Appendix

Bluephase G4 is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 3th (see the following tables) and IEC 60601-1-2:2014 4th Edition (according clause 7 and 8.9, tables 4 to 9). Bluephase G4 is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

## Warning

Portable wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie etc. can affect the Bluephase G4 and should be kept at least a distance of 30 cm (12 inches) to any part of Bluephase G4.

## Guidance and manufacturer's declaration – electromagnetic emission

The following tables are guidelines according to the medical standard IEC 60601-1-2.

Bluephase G4 is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase G4 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Bluephase G4 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Bluephase G4 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	N/A

Table: According to IEC 60601-1-2

## Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase G4 is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase G4 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of Bluephase G4 requires continued operation during power mains interruptions, it is recommended that Bluephase G4 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

Table: According to IEC60601-1-2

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase G4 is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase G4 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of Bluephase G4, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			<b>Recommended separation distance</b>
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

Table: According to IEC60601-1-2

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Bluephase G4 is used exceeds the applicable RF compliance level above, Bluephase G4 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Bluephase G4.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and "Bluephase Style G4"

Bluephase G4 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Bluephase G4 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Bluephase G4 as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.12	0.23
0.1	0.11	0.38	0.73
1	0.35	1.2	2.3
10	1.1	3.8	7.3
100	3.5	12	23

Table: According to IEC60601-1-2

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Ivoclar Vivadent – worldwide

**Ivoclar Vivadent AG**  
Benderstrasse 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Pty. Ltd.**  
1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 9795 9599  
Fax +61 3 9795 9645  
[www.ivoclarvivadent.com.au](http://www.ivoclarvivadent.com.au)

**Ivoclar Vivadent GmbH**  
Tech Gate Vienna  
Donau-City-Strasse 1  
1220 Wien  
Austria  
Tel. +43 1 263 191 10  
Fax +43 1 263 191 111  
[www.ivoclarvivadent.at](http://www.ivoclarvivadent.at)

**Ivoclar Vivadent Ltda.**  
Alameda Caliapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
[www.ivoclarvivadent.com.br](http://www.ivoclarvivadent.com.br)

**Ivoclar Vivadent Inc.**  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Tel. +1 905 670 8499  
Fax +1 905 670 3102  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

**Ivoclar Vivadent Shanghai Trading Co., Ltd.**  
2/F Building 1, 881 Wuding Road,  
Jing An District  
200040 Shanghai  
China  
Tel. +86 21 6032 1657  
Fax +86 21 6176 0968  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Marketing Ltd.**  
Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 3399  
Fax +57 1 633 1663  
[www.ivoclarvivadent.co](http://www.ivoclarvivadent.co)

**Ivoclar Vivadent SAS**  
B.P. 118  
74410 Saint-Jorioz  
France  
Tel. +33 4 50 88 64 00  
Fax +33 4 50 68 91 52  
[www.ivoclarvivadent.fr](http://www.ivoclarvivadent.fr)

**Ivoclar Vivadent GmbH**  
Dr. Adolf-Schneider-Str. 2  
73479 Ellwangen, Jagst  
Germany  
Tel. +49 7961 889 0  
Fax +49 7961 6326  
[www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

**Ivoclar Vivadent Marketing (India) Pvt. Ltd.**  
503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 22 2673 0302  
Fax +91 22 2673 0301  
[www.ivoclarvivadent.in](http://www.ivoclarvivadent.in)

**Ivoclar Vivadent Marketing Ltd.**  
The Icon  
Horizon Broadway BSD  
Block M5 No. 1  
Kecamatan Cisauk Kelurahan  
Sampora  
15345 Tangerang Selatan – Banten  
Indonesia  
Tel. +62 21 3003 2932  
Fax +62 21 3003 2934  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent s.r.l.**  
Via del Lavoro 47  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 6113555  
Fax +39 051 6113565  
[www.ivoclarvivadent.it](http://www.ivoclarvivadent.it)

**Ivoclar Vivadent K.K.**  
1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
[www.ivoclarvivadent.jp](http://www.ivoclarvivadent.jp)

**Ivoclar Vivadent Ltd.**  
4F TAMIYA Bldg.  
215 Baumoe-ro, Seocho-gu  
Seoul, 06740  
Republic of Korea  
Tel. +82 2 536 0714  
Fax +82 2 6499 0744  
[www.ivoclarvivadent.co.kr](http://www.ivoclarvivadent.co.kr)

**Ivoclar Vivadent S.A. de C.V.**  
Calzada de Tlalpan 564,  
Col Moderna, Del Benito Juárez  
03810 México, D.F.  
México  
Tel. +52 (55) 50 62 10 00  
Fax +52 (55) 50 62 10 29  
[www.ivoclarvivadent.com.mx](http://www.ivoclarvivadent.com.mx)

**Ivoclar Vivadent BV**  
De Fruittuinen 32  
2132 NZ Hoofddorp  
Netherlands  
Tel. +31 23 529 3791  
Fax +31 23 555 4504  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Ltd.**  
12 Omega St, Rosedale  
PO Box 303011 North Harbour  
Auckland 0751  
New Zealand  
Tel. +64 9 914 9999  
Fax +64 9 914 9990  
[www.ivoclarvivadent.co.nz](http://www.ivoclarvivadent.co.nz)

**Ivoclar Vivadent Polska Sp. z o.o.**  
Al. Jana Pawła II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 5496  
Fax +48 22 635 5469  
[www.ivoclarvivadent.pl](http://www.ivoclarvivadent.pl)

**Ivoclar Vivadent LLC**  
Prospekt Andropova 18 korp. 6/  
office 10-06, 115432 Moscow  
Russia  
Tel. +7 499 418 0300  
Fax +7 499 418 0310  
[www.ivoclarvivadent.ru](http://www.ivoclarvivadent.ru)

**Ivoclar Vivadent Marketing Ltd.**  
Qlaya Main St.  
Siricon Building No.14, 2<sup>nd</sup> Floor  
Office No. 204, P.O. Box 300146  
Riyadh 11372  
Saudi Arabia  
Tel. +966 11 293 8345  
Fax +966 11 293 8344  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent S.L.U.**  
Carretera de Fuencarral nº24  
Portal 1 – Planta Baja  
28108-Alcobendas (Madrid)  
Spain  
Tel. +34 91 375 78 20  
Fax +34 91 375 78 38  
[www.ivoclarvivadent.es](http://www.ivoclarvivadent.es)

**Ivoclar Vivadent AB**  
Dalvägen 14  
169 56 Solna  
Sweden  
Tel. +46 8 514 939 30  
Fax +46 8 514 939 40  
[www.ivoclarvivadent.se](http://www.ivoclarvivadent.se)

**Ivoclar Vivadent Liaison Office**  
: Tesvikiye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 0802  
Fax +90 212 343 0842  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Limited**  
Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SD  
United Kingdom  
Tel. +44 116 284 7880  
Fax +44 116 284 7881  
[www.ivoclarvivadent.co.uk](http://www.ivoclarvivadent.co.uk)

**Ivoclar Vivadent, Inc.**  
175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

691468

Date information prepared: 2018-06 / Rev. 0

691468/WE1

**ivoclar**  
**vivadent**  
*clinical*