REF BA110200

9795065

Ultimate BASE290

LED Curing Light

Lámpara LED de fotocurado
Lampe à Polymériser LED
LED-Polymerisationslampe
Lampada polimerizzante a LED
Led-uthardingslamp

MODEL: BASE290

P/N: BA110200

INSTRUCTIONS FOR USE
MODO DE EMPLEO
MANUEL D'UTILISATION
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING

B.A. INTERNATIONAL LTD.
Unit 9, Kingsthorpe Business
Centre Studland Road,
Northampton - NN2 6NE,
UK

Rev. 2018/10

Disclaimer:

1. First of all, we would like to thank you for purchasing BA Ultimate BASE290 LED Curing Light. Please carefully read the user manual and operate the device accordingly. Please keep this user manual well preserved for future reference.
2. BA International shall not assume any responsibility for any malfunction or damage that was caused by improper operation, removal, modification or maintenance.
3. We shall not assume any responsibility for malfunction or damage that was caused by operating in an unsuitable environment.
4. We cannot be held liable for any malfunction or damage caused by improper use, poor maintenance or force majeure.
5. BA International reserves the right to decline any responsibility for defective products which may have been caused by natural disasters, abnormal temperature or humidity that is inconsistent with the instructions in this manual.

Content:

• Intended Use.....	3
• Explanation of Graphical Symbols.....	3
• Pack Contents.....	4
• Instructions for Use.....	4
• Getting Started.....	4
• Technical Information.....	8
• Warnings & Safety Notes.....	9

Intended Use:

BASE290 is a portable light source for polymerising photo activated resins. This unit is for use only in dental clinics and is to be operated only by trained personnel.

Explanation of Graphical Symbols:

Storage conditions: -20°C to 60°C

CE Logo

Unit contains battery and electronic devices. These must be disposed of separately into a waste stream specifically intended for the recycling of such items.

Attention - consult accompanying documents

Centre Positive

UV Light

Refer to the instruction manual.

Class II electrical equipment

Applied part classification - Type B.

Manufacturer

Manufacturing Date

Pack Contents:

1x Power supply
1x Multi-plug (UK/EU/AUS)
1x Charger stand
1x Cure Unit
1x Pack of 50 protection sheaths
1x Protection shield

Instruction for use:

For your safety and the safety of your patients it is essential that you read and understand these instructions for use and safety notes.

Getting Started:

- Select the multi-plug for your location and fit to the power supply.
- Plug the power supply into a 100-240VAC outlet.
- Plug the DC plug into the charger stand socket.

CHARGING

Place the BASE290 cure unit into the charger stand.


A rotating display indicates charging is in progress - the light string length increases with charge level.

Constant illuminated display with an occasional dip in intensity indicates full charge.

Full charge should be reached in approximately four hours.


Note: The unit cannot be over charged.

NAVIKEY INTERFACE


Activating the light, selecting the mode of operation and user settings are all controlled via the Navikey.

The LED indicators adjacent to the Navikey display the current operating mode.

Power settings are selected by moving the Navikey to the right or left and are displayed via the line of blue indicators.

Duration of cure settings are selected by moving the Naviskey forwards and backwards are displayed via the row of green indicators.

The cure cycle is activated by pressing the Naviskey.

A second press of the Naviskey will terminate the cure cycle at any time.

On power levels three and four, a three second delay follows the cure cycle with an alarm to indicate that re-activation is available.

End of cure cycle is indicated by an audible and/or tactile (vibrate) alarm, both of which are optional.

To enable/disable the audible alarm; set the cure power to minimum setting and then press and hold the key to the left for three seconds, the unit will confirm a change with a bleep.

To enable/disable the tactile alarm; set the cure power to maximum setting and then press and hold the key to the right for three seconds, the unit will confirm a change with a vibrate.

If a cure cycle is interrupted by all lights flashing and the unit vibrating 7 times, this indicates that the battery level has dropped too low for the handpiece to operate.

MODES OF OPERATION

HIGH POWER

Power Setting	Power (P) (mW/cm ²)	Duration (T) (Seconds)			
		1	3	5	
Level 4	3000	1	3	5	
Level 3	2000	1	3	5	
Level 2	1000	1	5	10	20
Level 1	600	1	5	10	20

To assure correct curing follow the composite manufacturer's individual product guidelines for cure power and duration.

RAMP UP MODE

The RAMP UP mode is applicable to polymerize indirect resin composite restorations, light cure resin-based temporaries and other dental prosthesis or denture repair made in a dental laboratory.

BASE 290 cure unit incorporates a ramp up mode, starting at 600 mWcm⁻² and ramping up to 1000 mWcm⁻² over a 10 second period, followed by an alarm and a further 1000 mWcm⁻² to 2000 mWcm⁻² over a further 10 second period.

This function is selected by setting the power level to maximum followed by pressing and holding the Navkey. The power indicators will indicate ramp function is in operation by flashing.

WARNING:

As it's the case with all high performance curing lights, the high light intensity combined with long or repetitive curing cycle can generate heat. Prolonged exposure of areas near pulp and soft tissues could result in irreversible damages. Therefore, this product must only be operated by trained professionals clinicians.

POWER SAVING MODE

BASE 290 cure unit will revert to power saving mode after 2 minutes of inactivity. The previously set modes will be recalled following any press of the Navkey.

MAINTENANCE AND CLEANING

The single use BASE 290 barrier sheaths should be used and disposed of. The light unit should then be cleaned and disinfected using a cleaning cloth moistened with disinfectant.

Do not immerse the unit in water.

Follow local directives and guidelines for cleaning and disinfecting.

Use disinfectants which do not contain chlorine and are free from protein fixing effects e.g. chlorhexadine and always follow the manufacturer's guidelines.

Check the lens of the unit for cured dental resin; remove carefully paying attention not to scratch the lens.

WARNING: NOT AUTOCLAVABLE

Periodically check the power output with a light meter in the range 380nm to 500nm.

TECHNICAL INFORMATION

BASE290 CURE UNIT:

Classification	Class 1 medical device
Wavelength range	380nm-500nm
Wavelength peaks	395nm, 460nm
Light power	600,1000,2000 and 3000mWcm ⁻²
Weight	152g
Applied Part	Type B
Lens surface area	95mm ²

CHARGER STAND

Input	6V DC, 1.5A
Weight	587g

POWER SUPPLY

Input	100-240VAC 50/60Hz 0.5A
Max Output	6V DC, 1.5A

ENVIRONMENTAL CONDITIONS

Operation:	
Temperature	5°C to 30°C
Altitude	-390m to 3000m
Atmospheric Pressure	70kPa to 106kPa
Relative Humidity	15% to 90% non-condensing

TRANSPORT AND STORAGE:

Temperature	-20°C to 60°C
Altitude	-390m to 5'500m
Atmospheric Pressure	50kPa to 106kPa
Relative Humidity	8% to 80% non-condensing at 40°C

ACCESSORIES

BA110540 BASE 290 BARRIER SHEATHS(50pack)

BA110545 **BASE 290 PROTECTION SHIELD**

SERVICE

The BASE 290 cure unit is non field repairable.
If there is a problem please return to the distributor.

WARRANTY

As manufacturer, B.A. International is liable for material or manufacturing defects within a warranty period of one year from the date of purchase. B.A. International accepts no responsibility for damage caused by incorrect handling or by repairs carried out by third parties not authorised to do so by B.A. International.
Claims under warranty - accompanied by proof of purchase - must be sent to the vendor. The provision of service under warranty extends neither the warranty period nor any other guarantee period.

WARNINGS & SAFETY NOTES

- In accordance with EN60601-1, the unit is not suitable for use in potentially explosive atmosphere or with potentially explosive mixture of anaesthetic substances containing oxygen or nitrous oxide. This unit is for use in a dental clinical environment only.
- Electromagnetic interference: the unit has been tested and approved for all the current standards (EN 60601-1-2) both with regard to emission (electrical field rays and supply voltage interference) and immunity (electrical fields rays, intermittent bursts on the supply, electrostatic discharges and supply shock waves).
- Do not look directly at the light source and do not stare at light reflected from the tooth surface. Prolonged exposure may cause retinal damage.
- Restrict the exposure to the area of the aural cavity in which treatment is intended.
- Prevent heat accumulation when polymerizing with a rubber dam in place.
- Suitable blue light filtering safety glasses must be worn by patient and dentist during use.
- BASE 290 curing light emits an extremely high intensity light with an energy output very similar to high intensity quartz halogen lights. One should NOT direct light energy or place lens directly on soft tissues, oral mucous membrane, skin or the pulp for extended periods of time. (We recommend such for standard halogen lights as well.)
- Red or pink colored tissue can absorb more blue light energy than a white reflective tooth and/or tooth colored restoration. Limit contact with soft tissue as much as possible.
- Maintain a minimum of two millimeters between the tip of the light and the soft tissue. If gingival tissue is in close proximity to the area being cured i.e., Class II or Class V composite restorations, DO NOT expose gingival tissue for more than twenty seconds. A twenty-second cure is adequate (even when directing light interproximally for Class II restoration).
- Any unit used for polymerization, including the BASE 290 cure unit, must not be used on, or by, any person who suffers, or has suffered in the past, from photo-biological reactions (including urticaria solaris or protoporphyrie erythropetique), or any person receiving medical treatment and using photosensitive medication (including methoxsalenes or dymethylchlorotetracycline).
- Any person who has previously suffered from a retina or crystalline condition or who has undergone eye surgery, in particular cataract treatment, should visit an eye specialist before using the unit or before undergoing treatment with the unit.
- Under normal conditions, the power of the LED source does not vary and consequently the power level does not need to be checked as with ordinary curing units. However, if in doubt, use a light meter.
- Condensation inside an electrical appliance can be dangerous. If the unit needs to be transported from a cold to a warm area, it should not be used immediately but should be left long enough to reach the new room temperature.
- Do not use the unit near to any heat source. Use of solvent or flammable products on this unit could cause damage.
- For use only by a dental professional within a clinical environment.
- To assure correct curing, follow the composite manufacturer's individual product guidelines for cure power and duration.

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