INSTRUCTIONS FOR USE

Product description:

Etching gel is a 37% phosphoric acid etching gel with optimum thixotropic flow properties. It is used to etch tooth substances during the enamel-etch and total-etch techniques. If the safety instructions and directions on use are taken into account the product has no side-effects, according to the current state of knowledge. Underfilling materials can be etched by the gel in the event of longer application times.

- Only for dental purposes.
- Do not use after the end of the expiry date.
- Close the refill syringe immediately after use.
- Application syringes and application cannulas may only be used once.

Directions for use:

Filling the empty syringe applicators:

An empty syringe applicator is attached to the supply syringe with an adaptor (turn clockwise 180°). The small syringe is filled by carefully and slowly pressing the large syringe plunger. The small syringe plunger may also be pulled at the same time. Make sure that the small syringe is not overfilled or the plunger pushed out of the syringe! Then the small syringe applicator is unscrewed from the adaptor and the 60g syringe is closed and sealed again. A cannula is applied to the syringe applicator.

The etching gel can now be conveniently applied by carefully pressing the syringe stamp. The filled application syringe and the application cannula used are for single patient use only, the syringe must not be refilled. When using accessory application syringes, the expiry date and the LOT number of the storage syringe must be manually transferred to the label of the filled application syringe.

Preparation:

The edges of the areas of the tooth substance to be etched have to be clearly marked.

1.) The use of a dental dam is recommended! Any contact with soft tissue (skin, oral mucosa, eyes etc.) is to be avoided. If etching gel touches soft tissue, rinse the area immediately with plenty of water and where necessary contact a doctor.

2.) Before preparation, clean the tooth using a rubber cup with fluoride and fat-free polishing paste.

3.) Carefully remove cleaning residue with water, clean the surface of the tooth with oil-free pressurised air and then dry.

4.) Prepare the cavity properly, bevel the enamel margins to achieve permanent filling edges. Rinse the cavity well with water and dry with oil-free pressurised air. Avoid overdrying!

5.) Protect the dental pulp in the deep cavities (e.g. cover with a thin layer of a suitable underfilling material). The cavity should be dry and free of contamination while the overdrying of the dentine should be avoided.

Etching process:

6.) Apply the etching gel on the enamel edges of the cavity and, where necessary, to the dentine.

Recommended application time:

15 - 30 sec. on enamel

10 - 15 sec. on dentine

50 - 80 sec. on milk teeth and teeth with high fluoride content

7.) Rinse for at least 20 sec. with water.

8.) Dry with oil and water-free compressed air. The etched surface is not allowed to be touched until the treatment is continued, or be contaminated by saliva. If, however, there is contamination after etching and drying, the etching has to be repeated (approx. 10 sec.).

The etched enamel should look chalky white!

After the etching you can continue with the application of adhesives according to the manufacturer's directions for use.

Safety instruction:

H 314: Causes severe skin burns and eye damage. **Contains: Orthophosphoric acid.** Wear protective gloves/eye protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN: Gently wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor. Store in a closed container. Dispose of contents/container to

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According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

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