# EuroSept® Xtra Surface Disinfection System Wipes M

- (EN) Disinfection wipes
- (ES) Toallitas desinfectantes
- (FR) Lingettes de désinfection
- (DE) HS-Desinfektionstücher
- (IT) Salviette disinfettanti
- (NL) Desinfectieoekjes
- (PT) Pano de desinfeção
- (CS) Dezinfekční ubrousky
- (PL) Chusteczki dezynfekcyjne
- (DA) Desinficeringsklude
- (SV) Desinfektionsservetter
- (NO) Desinfiserende våtservietter

see imprint / ver impresión / voir l'impression / siehe Aufdruck / vedere l'impronta / zie afdruk / ver impressão viz otisk / patrz nadruk / se aftryk/ se avtrykk/ se avtrykk

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
EN: **Type of application:** As a medical device: Ready-to-use wet wipes for prophylactic surface disinfection of non-invasive medical devices in the medical field. **Notes on application:** Peel the stick & lift label from the refill bag and fix it to the dispenser bucket. The opening and shelf life dates must be written manually in the space provided. Open the bag top of the bag and places it in the dispensing bucket. Pull the first wipe of the fleece roll and thread it through the star shaped opening of the lid. Press the lid firmly back onto the container and close tightly. Carefully wipe the surface to be disinfected and leave to take effect. All surfaces must be completely covered to ensure disinfection. Always close the container after. After opening the refill bag and placing it in the dispenser, the shelf-life of the wipes is 28 days. The reprocessing of the dispensing bucket is to be carried out according to the adjacent reprocessing instructions using a suitable disinfectant wipe. Very good material compatibility, even with materials that are otherwise incompatible with alcohol. Suitable for the medical sector (dental practices, hospitals, nursing homes, care and rescue services) and the food, industrial, household and public sector. **Application times of the liquid:** Bacteria as per EN 13727, EN 1276, EN 13697 (dirty conditions): Conc. 1 min. "Limited virucidal" (enveloped viruses such as HBV, HIV, HCV as per EN 14476, dirty conditions): Conc. 1 min. "Unlimited virucidal P1" (enveloped viruses such as HBV, HIV, HCV and noro-, ade- and rotaviruses as per EN 14476 (clean conditions): Conc. 1 min. Polymyxoniums (SV40) (RRV/DV with contamination): Conc. 1 min. **Application times of the impregnated wipes:** Bacteria and Candida as per EN 13697 (dirty conditions): Conc. 1 min. Limited virucidal P1 as per EN 14476 (dirty conditions): Conc. 1 min. 100g impregnated with EuroSept Xtra Surface Disinfectant. Safety instructions: 228 Flammable solid. H313 Causes serious eye irritation. H336 May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/ attention. Store in a well-ventilated place. Keep container tightly closed. Dispose of contents/ container in accordance with all local, regional, national and international regulations. **100g impregnation liquid contain:** 20.0g 1-propanol, 5g 2-propanol. Recycle empty packaging only. Dispose of any remaining filling material as per EAK 07 06 99. (EC) 648/2004 (impregnation liquid): Contents: fragrances. All serious incidents related to the device shall be reported by the user/ patient to the manufacturer and to the competent authority of the State in which the user/ patient is established.

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PRISMAN GmbH