1 Intended use**Intended purpose**

Bonding of dental restorations to dental hard tissue, conditioning of Tetric CAD restorations, sealing of dentin lesions

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required

Application

For dental use only.

Description

Adhese® Universal is a light-curing, single-component dental adhesive for enamel and dentin. It is compatible with all etching techniques (self-etch, selective-enamel-etch and etch & rinse techniques). Adhese Universal can either be light-cured for 10 s at a light intensity of 500 to 1,400 mW/cm² or for 5 s at a light intensity of 1,800 to 2,200 mW/cm². In the restoration of Class I and II cavities, Adhese Universal can be light-cured in only 3 s from the occlusal aspect at a light intensity of 2,700 to 3,300 mW/cm² (e.g. using Bluephase PowerCure).

Indications

Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region

Contraindication

Do not use Adhese Universal if the patient is known to be allergic to any of the material's ingredients.

Limitations of use

- Adhese Universal must not be used for direct pulp cappings.
- Do not use Adhese Universal as a primer for ceramic restorative materials. A suitable ceramic primer must be used (e.g. Monobond Etch & Prime®).
- Do not use Adhese Universal for applications in which sufficient illumination cannot be ensured (e.g. cementation of root canal posts).
- Adhese Universal must not be used if the stipulated working technique cannot be employed.

Side effects

In rare cases, components of Adhese Universal may lead to sensitization. In such cases, Adhese Universal should no longer be used.

Interactions

- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used ≤ 14 days prior to the application of Adhese Universal.

Adhese® Universal

EN Instructions for Use Light-curing single-component dental adhesive	FI Käyttöohjeet Valokovettainen yksikomponenttinen sidosaaine
DE Gebrauchsinformation Lichthärtendes Einkomponenten-Dentaladhäsiv	NO Bruksanvisning Lyshefende enkomponent-dentaladhesiv
FR Mode d'emploi Adhésif dentaire monocomposant photopolymérisable	NL Gebruiksaanwijzing Licht uithardend ééncomponentadhesief voor tandheelkundig gebruik
IT Istruzioni d'uso Adesivo dentale monocomponente	EL Οδηγίες Χρήσεως Φωτοπολυμερίζόμενο οδοντιατρικό συγκολλητικό ενός φασκιδίου
ES Instrucciones de uso Adhesivo dental mono-componente y fotopolimerizable	TR Kullanma Talimatı Işıkla sertleşen tek bileşenli dental adeziv
PT Instruções de Uso Adesivo dental fotopolimerizável monocomponente	RU Инструкция по применению Светоотверждаемый однокомпонентный стоматологический адгезив

Phenolic substances (e.g. eugenol) may inhibit the polymerization of Adhese Universal.

- Before any treatment with Adhese Universal, completely rinse off all iron- and aluminium-based styptic agents with water.
- Fluoridation agents must not be used prior to the application of Adhese Universal.

Clinical benefit

- Reconstruction of chewing function
- Improvement of speech function
- Restoration of esthetics

Composition

HEMA, Bis-GMA, D3MA, MDP, ethanol, water, methacrylate-modified polyacrylic acid, silicon dioxide, camphor quinone, ethyl p-dimethyl aminobenzoate, 2-dimethyl aminoethyl methacrylate.

2 Application

2.1 Procedure for indirect and direct restorations

2.1.1 Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam Plus, is required.

2.1.2 Pulp protection / cavity liner

The 3sCure curing mode must not be used for restorations in areas close to the pulp. In very deep cavities, areas close to the pulp must be selectively coated with a calcium hydroxide liner (e.g. ApexCal®) and subsequently covered with a pressure-resistant cement (e.g. a glass ionomer cement).

2.1.3 Conditioning with phosphoric acid gel

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must always be conditioned with phosphoric acid gel (e.g. Total Etch). Please observe the instructions for use for the phosphoric acid gel.

a. Selective enamel etching

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15–30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

b. Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then onto the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

2.1.4 Handling the VivaPen® and bottle

a. VivaPen®

- Remove the cap of the VivaPen by turning it beyond the point where resistance is felt and simultaneously pulling it off.
- Attach the VivaPen brush cannula by snapping it into place. The attached cannula can be rotated to the desired position up to the point where resistance is felt.
Caution: Rotating the cannula beyond the point of resistance will detach the cannula!
- For hygienic reasons, always use a VivaPen protective sleeve by sliding the Pen with the cannula into the protective sleeve until the cannula punctures through the sleeve.
- Remove the paper separation strip before using the protective sleeve.
- The cannula tube can be bent up to 90°.
Caution: Bending the cannula several times or bending it at the cannula hub may break the cannula. Excessive bending may result in blocking of the cannula tube.
- A new cannula and protective sleeve must be used for each patient.
- Press the push button firmly down and wait until the cannula tip is saturated with Adhese Universal. Repeat if required. The cannula tip is saturated when the yellow colour of Adhese Universal is visible.
- If it is used several times or after a certain time on the same patient, the moistened cannula must be protected from light.
- Remove the protective sleeve and cannula after having used the VivaPen. The VivaPen brush cannula is removed by turning it beyond the point of resistance and simultaneously pulling it off.
- Close the VivaPen with the original cap after use.

Application notes

- When using the VivaPen for the first time, several clicks may be required to wet the cannula due to the engineering of the VivaPen. This does not influence the indicated fill-level or the function of the VivaPen.
- Do not immerse the VivaPen in disinfectant solutions or autoclave it.
- Contamination on the VivaPen can be removed, e.g. by wiping with alcohol or other non-oxidizing disinfectants.
- If the cannula has been contaminated (e.g. the cannula has dropped on the floor), it must be replaced by a new one.
- Cannulas of other manufacturers are not compatible with the VivaPen and therefore must not be used.
- The coloured indicator of the transparent fill-level gauge identifies the remaining amount of adhesive in the VivaPen. If the coloured indicator is no longer displayed, the VivaPen is empty.

b. Bottle

- Dispense the desired amount of Adhese Universal into a mixing well and apply it using a disposable applicator (e.g. Vivadent® Applicator Brush).
- Protect the adhesive from light (e.g. VivaPad®).
- Carefully close the bottle immediately after each use.

- A new disposable applicator must be used for each application.

Application notes

- Contamination on the bottle can be removed, e.g. by wiping with alcohol or other non-oxidizing disinfectants.
- Do not use the bottle intraorally. For hygienic reasons, suitable application aids should be used (e.g. VivaPad®).

2.1.5 Application of the adhesive

- Starting with the enamel, completely coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.

Important information:

- Avoid pooling, since this can compromise the fitting accuracy of the permanent restoration.
- Exposure to intensive light should be avoided during application.
- Prevent any contamination of the applicator/cannula with blood, saliva, sulcus fluid or water during application. In case of contamination, the application must be repeated with a new applicator/new cannula.

2.1.6 Light-curing the adhesive

Adhese Universal can either be light-cured for 10 s at a light intensity of 500 to 1,400 mW/cm² or for 5 s at a light intensity of 1,800 to 2,200 mW/cm². In the fabrication of Class I and II restorations, Adhese Universal can be light-cured in only 3 s from the occlusal aspect at a light intensity of 2,700 to 3,300 mW/cm² (e.g. using Bluephase PowerCure).

The instructions for use of the curing light must be observed.

Light intensity mW/cm ²	Exposure time
500 – 900	10 s
900 – 1400	10 s
1800 – 2200	5 s
2700 – 3300	3 s

Safety note

- Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.
- The 3sCure mode must not be used in case of caries profunda and very deep cavities.

2.1.7 Application of the restorative or luting composite

- Proceed according to the instructions for use of the restorative and/or luting composite used.

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ivoclar
vivadent®

2.2 Procedure for intraoral repairs of fractured composite fillings/ composite restorations with light-curing materials

- Roughen the surface to be repaired using diamond finishers and thoroughly clean with water spray and dry with oil- and water-free compressed air.
- Proceed according to the direct application for Adhese Universal.

2.3 Procedure for the conditioning of Tetric® CAD restorations

- Proceed according to the instructions for use for Tetric CAD.

3 Safety information

- Adhese Universal is irritant. Avoid contact with the skin, mucous membrane and eyes.
- After skin contact, immediately rinse with copious amounts of water.
- If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ ophthalmologist.
- In rare cases, contact with the skin may lead to sensitization towards the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).
- The current Summary of Safety and Clinical Performance (SSCP) can be found in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

Warnings for phosphoric acid gel:

Phosphoric acid gel is corrosive. Avoid contact with the eyes, mucous membrane and skin (protective goggles are recommended for both patients and operators). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Failure of the adhesive bond (loss of the restoration)
- Postoperative sensitivity

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on packages

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.