REF BASE 800

BA OPTIMA U-MIX

Capsule Mixer

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

UNIVERSAL KAPSELMISCHER

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MADE IN TAIWAN
Fabricado en Taiwan
Fabriqué en Taiwan
Hergestellt in Taiwan

B.A. INTERNATIONAL LTD.
Unit 9, Kingsthorpe Business Centre
Studland Road, Northampton - NN2 6NE
U.K.

Rev. 2018/07

BA OPTIMA U-MIX

Disclaimer:

1. First of all, we would like to thank you for purchasing BA OPTIMA U-MIX Capsule Mixer. Please carefully read the user manual and operate the apparatus accordingly. Please keep this user manual well preserved for future reference.
2. BA shall not assume any responsibility for any malfunction or damage that was caused by improper operation, removal, modification or maintenance.
3. We shall not assume any responsibility for malfunction or damage that was caused by operating in an unsuitable environment.
4. We cannot be held liable for any malfunction or damage caused by improper use, poor maintenance or force majeure
5. BA reserves the right to decline any responsibility for defective products which may have been caused by natural disasters, abnormal temperature or humidity that is inconsistent with the instructions in this manual.

Content:

• Intended Use.....	3
• Professional Use Statement.....	3
• Safety Precautions.....	4
• Explanation of Graphical Symbols.....	5
• Product parts & accessories.....	6
• Important tips before operation.....	7
• Operating Steps.....	8
• Important tips on product maintenance.....	9
• Important tips on product cleaning.....	9
• Troubleshooting.....	10
• Product Specification.....	11
• EMC Tables.....	12

Intended use:

"BA OPTIMA U-MIX Capsule Mixer" is intended for the mixing of dental material capsules. Applies to all capsules products on the market.

Professional Use Statement:

This device is capsule mixer for dental use. Used only in dental clinics and operated by trained personnel.

SAFETY PRECAUTIONS

Please follow the instructions in this manual for correct and safe operation of the device. Pay special attention to the following warning signs beside the operation descriptions where applicable.

DANGER :

This message appears where improper operation may cause severe injury to the human body if instructions are not followed correctly.

WARNING :

This message appears where improper operation may cause serious damage or defect to other objects or the device itself if instructions are not followed correctly.

CAUTION :

This message appears where improper operation may cause slight injury to the human body if instructions are not followed correctly.

IMPORTANT :

This message appears where improper operation may cause slight damage or defect to other objects or the device itself if instructions are not followed correctly.

EU-wide legislation, as implemented in each Member State, requires that waste electrical and electronic products carrying the mark (left) must be disposed of separately from normal household waste. This includes monitors and electrical accessories, such as signal cables or power cords. When you need to dispose of your display products, please follow the guidance of your local authority, or ask the shop where you purchased the product, or if applicable, follow any agreements made between yourself.

The mark on electrical and electronic products only applies to the current European Union Member States.

Explanation of Graphical Symbols:

 IEC 60878 safety Sign ISO 7010-M002:Refer to the instruction manual.

 Indicates the polarity of the DC power input.

 CLASS I equipment has met requirements of European Medical Device Directive 93/42/EEC.

Product Cleaning:

If there is any dirt or stains on the surface, please use a cloth to wipe clean. If it is hard to clean, please use a neutral cleansing agent. 70% isopropyl alcohol is recommended.

 **IMPORTANT :**
Never use a solvent, volatile oil or immerse for cleaning.

 **WARNING :**
Do not modify this device without authorization of the manufacturer.


 **IMPORTANT :**
Do not use other adapters without authorization from the manufacturer.

 **CAUTION :**
Please check the specification before using the adapter, the incorrectness could cause the malfunction of the product.

 **CAUTION :**
Please use dental material capsules. Do not use corrosive materials to cause the malfunction of the product.

Product Parts & Accessories

Please read the following information thoroughly before first use.


1. Protective Shield:
The safety switch will be activated when the Protective Shield is not closed. The machine will automatically stop its operation in case the protective shield is opened.
2. Holding Arms:
Where to place the capsule.
3. Power Button:
The Mode indicator & Time Display indicator will start to blink (showing display of "10") when the power button is pressed ON.
4. Mode Set-up Buttons:
Press Mode button to set required Mode setting.
5. Mode Indicator:
When the required mode is selected, the corresponding indicator will be turned on.
6. ▲ ▼ Timer Setting Button:
To set the mixing time, adjustable from 1-30 seconds.
7. Time Display indicator:
Displays the selected mixing time. The timer will blink during operation as it counts down.
8. START/STOP Button:
Press the button to start operation. Press once again to stop the operation at anytime.
9. Adapter(UK/EU/AUS) : 100V ~ 240V 50/60Hz 2A

Important tips before Operation


Due to transportation requirements, all the motors inside each unit are secured with 2 fixing screw bolts. Please observe the diagram below for removing these screws before placing the capsule mixer on the cabinet.

1. Please unscrew the two screw bolts, which are pre-installed at the bottom of each unit, as shown in the Diagram 1.


Please maintain the balance of the apparatus while removing the screw bolts. Please remove the two bolts simultaneously as unscrewing 1 bolt at a time could result in difficulties loosening the second bolt.

2. Please retain the two bolts to the screw holes as shown in Diagram 2, for future needs. When the apparatus is required to be transported from one place to another, please screw back these two screw bolts before removal or delivery to prevent possible damage.

3. Insert the "L" shape plug into the amalgamator power socket, then connect the power plug with the adaptor. Please make sure the power plug is connected securely (As shown in Diagram 3)

Operating STEPS

1. After pressing the Power button on, the mode indicator and the time indicator will start blinking. As shown in diagram 4.

2. Depending on each capsules different characteristics, the time required for mixing may vary from each other.
 - ※ The Amalgam capsule & Filling material require an average time of 10 seconds to gain a consistent mixture. Material with special specifications may vary with different capsules.
 - ※ Mixing time can be reset at any given time setting between 1 to 30 seconds. The user is able to setup two different time settings using the memory. Regarding individual capsule time setting, please review the instruction manual accompanied with the capsule products for further details.
- 1) Please make sure the Protective glass is securely closed.
 - 2) Press Set button to select the designated mode.
 - 3) Press ▲/▼ to change the time set up mode. Press ▲ or “▼” to add or reduce the mixing time by one second.
 - 4) Press “SET” button to confirm the setting. The time indicator will change from a blinking display to a fixed display.
3. Open the Protective shield, and place the capsule into the holding arms as shown in diagram 5.

4. Securely close the protective glass and press “START/STOP” button to begin mixing operation. The timer display on the screen will start to countdown as time proceeds.
 - This apparatus will not operate if the protective glass is not securely shut. The time indicator will display “OFF” as a result.
5. After Mixing is complete, please open the protective glass and take the capsule from the holding arms. Please review diagram 6 for removing the capsule.

Important tips on Product maintenance

1. Do not operate when the capsule is not in place. Operation without the capsule may result in malfunction. This damage is not covered under the warranty.
2. This product should not be placed in direct sunlight. DO NOT place this product in a dusty environment.
3. If any liquid escapes from the capsule to the inside of the mixing chamber, immediately wipe clean to prevent possible damage, colour change or metallic distortion caused by dental use liquid or solvent.

Important tips on Product cleaning

Clean the surface of this equipment with a soft clean cloth soaked in warm water or a mild neutral cleaning agent. Avoid using a volatile solvent when cleaning.

Trouble Shooting

In the unlikely event of a problem occurring with the capsule mixer, please check against the below listed SOP. If the solution is not provided below, please contact your local dealer or BA International for a solution or repair. Thank you.

Scenario	Possible Cause	Possible Solution
Unable to operate when pressing START / STOP button	The power adaptor might be loose.	Ensure the power adaptor is plugged in to the unit.
	The L shaped plower plug might be loose.	Ensure the L shaped plug is properly plugged in.
	The protective shield is open.	Make sure the protective shield is firmly shut.
The apparatus moves during operation.	Please check if the apparatus is placed on uneven surface.	Please place a securely stable plastic pad under the equipment.
Abnormal Noise	2 Screw bolts have not been removed.	Unscrew the two screw bolts at the bottom of the units.

Product Warranty & Repair

All B.A. International Products are manufactured under thorough quality control. If the apparatus were to be broken under normal operation in accordance with the Instruction Manual and within the warranty period, the broken unit is subject to repair and maintenance free of charge. Please contact your nearest local dealer or B.A. International for your inquiry & best solution. Our warranty period extends to 1 year from the date of purchase.

Product Specification

Product Classification:

Protection against power strike	CLASS I equipment
Voltage	DC24V
Input Power	50VA
Specification on the power adaptor	(Input) AC100V~240V 50/60Hz, 2A (Output) DC24V , 2.7A

System operational conditions	<p>Temperature: 25°C to 35°C (77°F~95°F)</p> <p>Humidity: 15% to 95% @ 40°C, noncondensing</p> <p>Atmospheric Pressure: 700-1060 hpa</p>
System storage conditions	<p>Temperature: -20°C to 65°C (-4°F~149°F)</p> <p>Humidity: 15% to 95% @ 40°C, noncondensing</p> <p>Atmospheric Pressure: 700-1060 hpa</p>
System transportation conditions	<p>Temperature: -20°C to 65°C (-4°F~149°F)</p> <p>Humidity: 15% to 95% @ 40°C, noncondensing</p> <p>Atmospheric Pressure: 700-1060 hpa</p>

EMC TABLES

Guidance and manufacturer's declaration – electromagnetic emissions		
The model BASE800 is intended for use in the electromagnetic environment specified below. The customer or the user of the model BASE800 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The model BASE800 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model BASE800 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Recommended separation distances between portable and mobile RF communications equipment and the model BASE800

The model BASE800 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model BASE800 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model BASE800 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The model BASE800 is intended for use in the electromagnetic environment specified below. The customer or the user of the model BASE800 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model BASE800 requires continued operation during power mains interruptions, it is recommended that the model BASE800 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE *UT* is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The model BASE800 is intended for use in the electromagnetic environment specified below.
The customer or the user of the model BASE800 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BASE800, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$ <p>.....</p> $d = 1,2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ <p>.....</p> $d = 2,3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 A/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model BASE800 is used exceeds the applicable RF compliance level above, the model BASE800 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model BASE800.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than V/m.

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EMC TABLES

Guidance and manufacturer's declaration – electromagnetic emissions

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Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The model BASE800 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model BASE800 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

**Recommended separation distances between
portable and mobile RF communications equipment and the model BASE800**

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BASE800, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 A/m	$d = 1,2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BASE800, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$ <p>.....</p> $d = 1,2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ <p>.....</p> $d = 2,3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 A/m	

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Recommended separation distances between portable and mobile RF communications equipment and the model BASE800			
The model BASE800 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model BASE800 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model BASE800 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The model BASE800 is intended for use in the electromagnetic environment specified below. The customer or the user of the model BASE800 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model BASE800 requires continued operation during power mains interruptions, it is recommended that the model BASE800 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE <i>UT</i> is the a.c. mains voltage prior to application of the test level.			

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BASE800, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 A/m	$d = 1,2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ <hr/> $d = 2,3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

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