

 HENRY SCHEIN®

 ZAHN™ LABORATORY
DIVISION

GIBRALTAR-ALABASTER

Plaster for Dental Cases, Type II

Escayola para Modelos, Tipo II

Plâtre pour Modèles, Type II

HS-Modellgips (Alabaster), Typ II

Gesso duro, Classe II



900-2878

White/Blanco/Blanc/Weiss/Bianco



For Professional Use Only
Unicamente Para Uso Profesional
Pour usage professionnel uniquement
Nur für den Dentalgebrauch
Solo per uso professionale

Contains:
Contiene:
Contient:
Inhalt :
Contenuto:

20kg

Technical data

Setting time/Tiempo de fraguado/Temps de prise/Gesamtaushärtezeit/Tempo di presa	10-12min
Compressive strength 1hr/Fuerza compresiva 1h/Resistance de compression 1 hr/Druckfestigkeit nach 1h/Resistenza alla compressione dopo 1 ora	15
Expansion 2 hrs/Expansión 2 horas/Expansion 2 heures/Expansion nach 2h/Espansioni di presa dopo 2 ore	0.15%
Mix ratio/Proporción de mezcla/Proportion de mélange/Mischungsverhältnis/Miscelatura	50ml/ 100g

MADE IN GERMANY

Fabricado en Alemania/Fabriqu  en Allemagne
In Deutschland hergestellt/Prodotto in Germania

Distributed by:
Distrib  par:

Distribuido por:
Vertrieb durch:
Distribuito da:

HENRY SCHEIN INC.
Melville, NY 11747 USA

Henry Schein U.K. Holdings Ltd.
Gillingham ME8 0SB U.K.