

BAGETTE ANWENDUNGSTABELLE

PART CODE	GRÖSSE	TÜTEN/BOX	PASSE INHALTE
IMS-1346	2-1/4" x 4" (57 mm x 102 mm)	200	Kofferdam-Klammern und Bohrer
IMS-1347	2-3/4" x 9" (70 mm x 229 mm)	200	1 Handstück oder 1-3 Handinstrumente
IMS-1236	3-1/2" x 9" (89 mm x 229 mm)	200	2-5 Handinstrumente
IMS-1348	5-1/4" x 10" (133 mm x 254 mm)	200	Handstück oder mehrere Instrumente
IMS-1237	7-1/2" x 13" (190 mm x 330 mm)	100	<ul style="list-style-type: none"> • 1/4 DIN Kasette – 5 Instrumente; IME14DIN5X • Exam Kasette – f. 5 Instrumente; IM605X • Double-Decker Kasette – f. 8 Instrumente; IM908X • Große Handstück Kasette; IM60HX • Kleine Ultraschallkasette; IM30HUX • Symmetry IQ® Handstück Kasette; IM60IQ8, IM30IQ8
IMS-1238	10" x 14" (254 mm x 356 mm)	100	<ul style="list-style-type: none"> • 1/2 DIN Kasette – f. 8-10 Instrumente; IME12DINXXX • Kleine Kasette – f. 8 Instrumente; IM508X • Kleine Kasette – f. 10 Instrumente; IM510X • Double-Decker Kasette – f. 14 Instrumente; IM914X • Radiographie Kasette; IM7000 • Space Saver Seven Kasette; IMSS7X • Kleine dünne KFO Kasette; IMORTHSM8 • Kleine Resin Kasette – f. 7-16 Instrumente; IMS-3107, IMS-3112, IMS-3116
IMS-1239	12" x 17" (305 mm x 432 mm)	100	<ul style="list-style-type: none"> • DIN Kasette – f. 13-20 Instrumente; IMEDINXXX • Große Kasette – f. 16 Instrumente; IM416X • Große Kasette – f. 20 Instrumente; IM420X • Double-Decker Kasette – f. 20 Instrumente; IM920X • Signa-Stat Kasette – f. 12 Instrumente; IM812X • Große dünne KFO Kasette; IMORTHLG8 • Große Resin Kasette – f. 12-24 Instrumente; IMS-1118, IMS-1124, IMS-1112, IMS-1116
IMS-1345	13" x 20" (330 mm x 508 mm)	100	<ul style="list-style-type: none"> • Chirurgie Kasette; IMOSX • Große Custom Kasette; 118

Besuchen Sie uns auf Hu-Friedy.eu

Customer Care Department: Free Phone 00800 48 37 43 39
Free Fax 00800 48 37 43 40 | Email: info@hu-friedy.eu

©2014 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-520HFE/1114

