

SOLDERING INVESTMENT

ES Soldadura

FR Revêtement à souder

DE HS-Lötmasse

IT Rivestimento per saldatura

NL Soldeer inbedmassa



Mixing ratio (powder : liquid)
 Proporción de mezcla (polvo: líquido)
 Rapport de mélange (poudre : liquide)
 Mischungsverhältnis (Pulver : Flüssigkeit)
 Rapporto di miscelazione (polvere: liquido)
 Mengverhouding (poeder : vloeistof)

100 g / 23 ml

Mixing time under vacuum
 El tiempo de mezclado al vacío
 Temps de malaxage sous vide
 Mischzeit unter Vakuum
 Tempo di miscelazione sotto vuoto
 Mengtijd onder vacuüm

30 sec

Mixing time
 El tiempo de mezclado
 Temps de malaxage
 Mischzeit
 Tempo di miscelazione
 Mengtijd

60 sec

Setting time
 El tiempo de curado
 Temps de durcissement
 Aushärtezeit
 Tempo di indurimento
 Uithardingtijd

Max 10 min

6 kg



+H65898809140U



Distributed by (in US only)

HENRY SCHEIN INC.
 135 DURYEA ROAD
 Melville, NY 11747 USA

EC REP

HENRY SCHEIN U.K. HOLDINGS LTD.
 MEDCARE HOUSE, CENTURION CLOSE
 GILLINGHAM BUSINESS PARK
 GILLINGHAM, ME8 0SB U.K.
www.henryscheinbrand.com

Rev. 2015/11