

Bluephase® G4

The curing light

Polyvision®
INSIDE 

With intelligent
curing assistant

Efficient
Esthetics


ivoclar
vivadent®
passion vision innovation

Intelligent LED light – carefree curing

Patented Polyvision technology is your personal assistant when it comes to curing dental materials.

This automatic assistance system helps you achieve reliable curing results. The curing light automatically detects movement of the handpiece during the light-curing procedure. It alerts the user of the improper operation by vibrating and automatically extends the exposure time by 10 per cent, if necessary. If the movement may prevent the material from curing properly, the light will automatically interrupt the exposure cycle so that it can be repeated. Bluephase G4 comes with active anti-glare protection to prevent the light from being activated in open space.

Additional advantages:

Reliable curing performance for

all light-curing materials due to polywave LED

Short curing times


starting from 10 seconds due to a high intensity of 1,200 mW/cm²

Wide 10-mm light guide for

time-saving single-exposure curing procedures

Additional Pre-Cure program

for the removal of excess cement



Contactless battery charging

Delivery forms

691258

Bluephase G4

691257

Bluephase G4 with Click & Cure and integrated radiometer

Homogeneous
light emission performance

Backlit display
for optimized hygiene

Light-curing assistant
Polyvision[®]

*The high-performance polymerization unit with
increased light intensity awaits you at IDS 2019.*

IDS
2019

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

702866/DE/2018-11



ivoclar
vivadent[®]
passion vision innovation